Schlauchpumpe
REGLO Analog

2 Kanäle
ISM830 6 Rollen
ISM829 8 Rollen
ISM795 12 Rollen

4 Kanäle
ISM828 6 Rollen
ISM827 8 Rollen
ISM796 12 Rollen

Tubing pump
REGLO Analog


2 channels
ISM830 6 rollers
ISM829 8 rollers
ISM795 12 rollers

4 channels
ISM828 6 rollers
ISM827 8 rollers
ISM796 12 rollers

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Please note

We recommend that you read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC® does not admit responsibility for the handling of chemicals.

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Safety precautions

ISMATEC® Pumps are designed for pumping applications in laboratories and industry. As such, it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed:

- ▶ The circuit between mains supply and pump has to be connected to earth ground.
- ▶ The pump must not be operated outside the designed operating and environmental conditions.
- ▶ The pump must not be used:
 - for medical applications on human beings
 - in explosion proof chambers or in the presence of flammable gases or fumes
- ▶ The pump must be switched OFF when pump-heads, cassettes or tubing are inserted or changed.

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Caution

Rotating pumphead creates a pinch and crush hazard.

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Safety precautions

- ▶ The permeability of tubing depends on the material used and pressure conditions. Tubing can also become electrostatically charged. Please be aware of possible hazards when routing tubing in explosion-proof chambers.
- ▶ Pump-heads consist of rotating parts. Therefore, the pump must not be operated before the cassettes are fully snapped-in.
- ▶ Do not manipulate the pump-head before the pump is switched OFF and disconnected from the mains supply.
- ▶ Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the rotating pump head.

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Safety precautions

- ▶ Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- ▶ Only new fuses, according to the specifications stated on Page 10 in this manual, must be used.
- ▶ The fuse-holder must not be short-circuited.
- ▶ Do not open or remove the housing while the pump is operating.
- ▶ Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- ▶ For service and repairs carried out by the customer or by third-party companies, ISMATEC® denies any responsibility.

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Warranty

For all parts manufactured by ISMATEC® from date of delivery: 2 years
All other parts, excluding consumables: 1 year

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

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Warranty terms

ISMATEC® warrants the perfect functioning of our products, provided they have been installed and operated correctly, according to our operating instructions. If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion.

A pump must be returned in the original ISMATEC® packaging or in packaging of equal quality. The duration of the warranty is not affected by making a claim for warranty service. Further claims are excluded. Shipping costs are charged to the customer.

Our warranty becomes invalid in the case of:

- ▶ improper operation by the user, or if the pump is diverted from its proper use
- ▶ unauthorized modification or misuse by the user or by a third party
- ▶ improper site preparation and maintenance
- ▶ operation outside of the environmental and electrical specifications for the product
- ▶ use of third-party software, hardware, accessories or consumables purchased by the user and which do not comply with our specifications.

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Product

The REGLO *Analog* tubing pump is available with 2 or 4 tube channels and 6, 8 or 12 pump rollers.

REGLO *Analog* 2 channels

(Speed 3.2–160 rpm)

Order Number	Type	Flow Rates ml/min
ISM830	MS-2/06	0.005–68
ISM829	MS-2/08	0.004–57
ISM795	MS-2/12	0.003–38

REGLO *Analog* 4 channels

(Speed 2.0–100 rpm)

Order Number	Type	Flow rates ml/min
ISM828	MS-4/06	0.003–43
ISM827	MS-4/08	0.003–35
ISM796	MS-4/12	0.002–24

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Product

Package contents

- ▶ Pump with 2 or 4 cassettes (type as ordered)
- ▶ Power Cord
- ▶ Operating Manual

Please check the packaging and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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Rear panel

1. Analog Interface
2. Mains Socket
3. Fuse-Holder with Voltage Selector 115/230V

Mains voltage

Mains voltage	Fuse rating
220–240 V _{AC} 110–120 V _{AC}	2 x T500mA/250V 2 x T500mA/250V

⚠ Socket/Power cord

Use exclusively the originally supplied power cord.

The socket must be connected to earth ground.

⚠ Before starting-up

Check that the voltage setting visible in the window of the fuse-holder complies with the local mains voltage. If necessary, the voltage setting must be changed and the 2 fuses must be replaced.

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
Voltage setting 115V / 230V and changing the fuses

➔ Switch the pump OFF, then disconnect the mains plug.

1. Remove the fuse-holder by opening it with a small screw-driver (size 0).
2. Remove the voltage selector plate. Turn it and re-insert it into the fuse-holder so that the required voltage rating is facing the window of the fuse-holder.
3. Insert 2 new fuses
230 V_{AC}: 2x T500mA/250V
115 V_{AC}: 2x T500mA/250V

⚠ Use always 2 slow-blow fuses of the same type complying with the local mains voltage.

4. Re-install the fuse-holder. The voltage rating is visible in the window.

⚠ Before starting-up

Check that the voltage setting visible in the window of the fuse-holder complies with the local mains voltage.

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


Bestell-Nr. IS3510A

Inserting the tubing

1. ➡ Switch the pump OFF
2. Remove the cassette by slightly pressing the fixing-tongue and lifting it simultaneously
3. Insert the 3-stop color-coded tubing with one collar into the cassette
4. Let the tubing hang down (prevent it from being twisted)
5. Insert the second collar of the tube at the other end of the cassette
6. Reinsert the cassette into the roller-head
7. ➡ Switch the pump ON

➡ **When the pump is idle, release all cassettes at the fixing-tongue (2).**

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Starting the pump

- ▶ Check the voltage rating in the window of the fuse-holder (on rear panel). If necessary, change the rating as indicated on Page 10–11.
- ▶ Connect the pump tubing to the system.
- ▶ Connect the pump to the mains.

 When using new tubing for the first time, it may occur that, depending on the tubing used (hardness and diameter), the pump cannot be primed and, hence, does not deliver the liquid. If that is the case it is recommended that the tubing be wet and run in the pump for 15 to 30 minutes.

➡ For applications with a high differential pressure (>1 bar), use of pressure lever cassettes are recommended.

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Starting the pump

1. Set the required speed on the 2-digit speed selector.
99 = 100 rpm or 160 rpm (max. revolution)
➔ The rotation speed can be adjusted while the pump is running.
2. Choose the rotation direction
3. Start the pump with the mains switch = run/stop function
4. Overload indicator (red LED)

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Overcurrent protector

The REGLO *Analog* tubing pump features an overload protector.

When an overload condition occurs, an activated overcurrent protector is indicated by the red LED (4) and the pump is stopped.

Whenever this situation occurs, **the pump must be switched OFF immediately.**

➡ Let the pump cool down (it takes about 2 minutes until the pump is ready again).

Before the pump is re-started, it is most important to check the reason for the overload (eg. too high differential pressure).

➡ Only when the cause of the overload has been detected and the failure corrected accordingly may the pump be started again.

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Pumping against pressure

The REGLO *Analog* can be used for continuous duty at a differential pressure of 1.0 bar max. for the smaller tube sizes. Pressure lever type cassettes allow operating at higher pressures.

In case of any doubts please contact your ISMATEC® agent.

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When the pump is not in use

When the pump is idle, pressure on the tubing should be released. Releasing the cassette on the right side is sufficient.

This helps to protect the tubing from unnecessary strain and prolongs its service-life.

Syphoning effect

When pressure on the tubing is released, the fluid can flow back to the reservoir.

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Running-in period for tubing

Every new tube requires a running-in period. If constant and reproducible flow rates are required, new tubing should be run-in with water or the medium to be pumped for 1 to 3 hours before the application is started.

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Please Note

Please refer to our detailed tubing documentation.

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Tubing life

The service-life of the tubing depends on the application and the tubing material used.

Example: chemicals, rotation speed, differential pressure, temperature, viscosity, pressure on tubing, etc.

General information on the service-life is stated in our detailed tubing documentation (without obligation!).

In order to improve the lubrication and service-life of the tubing, we recommend users to spray both the tubing and the pump rollers with our silicone oil spray (Order No. SC0179).

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
Digitale Eingänge (TTL-Pegel)

Digital inputs (TTL-level)

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Pin 2,	remote
Pin 3,	start
Pin 4,	direction
Pin 13,	internal speed

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Analog interface

Pin 1, GND (ground)

Reference potential for all other inputs.

Pin 2, remote

For changing between manual control and analog interface. For activating the analog interface, Pin 2 must be connected with Pin 1 (GND).

Pin 3, start

In remote operation (Pin 2 to GND) the pump starts when connected to Pin 1 (GND).

Pin 4, direction

In the **open** position, the pump turns counter-clockwise; when connected to Pin 1 (GND) it turns clockwise.

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
Analog interface

Pin 5, speed IN

External speed control
 (0–5V_{DC}, 0–10V_{DC}, 0–20mA, 4–20mA)
 Input impedance and input range can
 be selected via a dip-switch inside the
 pump (see P. 22).

Input impedance

0–5	V	18 kΩ
0–10	V	38 kΩ
0–20	mA	270 Ω
4–20	mA	270 Ω

Pin 7, +20V_{DC}

About +20 V_{DC} are available
 (max. current 0.2 A).

Pin 9, speed OUT

Frequency proportional to the rotation speed:
 2 channels: 0–8 kHz, 3.2–160 rpm
 4 channels: 0–5 kHz, 2.0–100 rpm

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
Analog interface

Pin 10, +5V_{DC}

About +5 V_{DC} are available (max. current 0.1 A).

Pin 13, internal speed

Analog interface activated (Pin 2 on GND)

- Pin 13 open: The rotation speed is adjusted via Pin 5 (speed IN).
- Pin 13 on GND: The rotation speed can be adjusted by the speed selector on the control panel of the pump.

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
Settings of switch S1

➔ Make sure that the pump is disconnected from the mains supply.

⚠ The instrument should only be opened by a qualified technician!

Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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Pins		Impedanz	DIP-Switch 1	DIP-Switch 2	DIP-Switch 3	DIP-Switch 4	DIP-Switch 5	DIP-Switch 6
Pin 5 speed IN	0-5V _{DC}	18kΩ	ON*	OFF*	OFF*	OFF*	OFF*	ON*
	0-10V _{DC}	38kΩ	OFF	ON	OFF	OFF	OFF	ON
	0-20mA	250Ω	OFF	OFF	ON	OFF	ON	ON
	4-20mA	250Ω	OFF	OFF	OFF	ON	ON	OFF

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Footswitch

Order No. ISM891

This footswitch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc. Both hands are free for handling the bottles and tubing.

3-stop tubing

For selecting the correct tubing, please refer to our detailed tubing documentation.

➔ When using new tubing for the first time, it may occur that, depending on the tubing used (hardness and diameter), the pump cannot be primed and, hence, does not deliver the liquid. If that is the case, the tubing should be primed and run for 15 to 30 minutes.

For applications with a high differential pressure (> 1 bar), pressure lever cassettes should be used.

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
MS/CA Anpresshebel
MS/CA pressure lever
MS/CA levier de pression

Zubehör

Ersatz-Kassetten

Kassette MS/CA Click'n'go Material
POM

Bestell-Nr. **IS3510A**

Kassette MS/CA Anpresshebel
Material **POM**

Bestell-Nr. **IS0649A**

Kassette MS/CA Anpresshebel
Material **PVDF**

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Accessories

Spare-cassettes

Cassette MS/CA Click'n'go Material
POM

Order No. **IS3510A**

Cassette MS/CA pressure lever
Material **POM**

Order No. **IS0649A**

Cassette MS/CA pressure lever
Material **PVDF**

Order No. **IS3629A**

➔ Pressure lever cassettes are available on request. This type of cassette may provide better results when using the Tygon MH tubing and at elevated differential pressure conditions.

➔ The cassettes allow the user to insert tubes with different diameters and materials on the same roller-head.

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*) Bestell-Nr. = Tygon® LMT-55

In the tables listed subsequently, the values indicated are only approximate and determined as follows: ml/min, per channel, with water and Tygon® tubing, without differential pressure.

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Please note

We will be pleased to help the user to select the tubing material. However, the user himself has the final responsibility for the selection of the correct tubing material.

For reproducible results we recommend you to use tubing from ISMATEC®.

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Disposal

Please retain packaging materials until the product warranty ends. Afterwards please discard packaging materials in an environmentally-friendly manner according to local regulations. Once the useful life of the product has ended, please ensure proper disposal according to local laws. Paper, cardboard, plastic, and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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Service and repairs

Like all pumps, the REGLO *Analog* contains components that will wear over a period of time.

For repairs, send the complete pump with a detailed description of the failure to your ISMATEC® distributor. Please use the original ISMATEC® packaging or packaging of equal quality.

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Technical Specifications

Drive

Motor type	DC-Motor
Speed range	
2 channels	3.2 – 160 rpm
4 channels	2.0 – 100 rpm adjustable in 1 % steps

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Please note

We also recommend you observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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Operating conditions

Temperature +5 to +40°C,
Rel. humidity max. 80%

– not condensing, at normal environmental conditions

Dimensions/Weight

(DxWxH)
2 channels 178x100x143 mm
(7"x4"x5 2/3")
4 channels 190x100x143 mm
(7 1/2"x4"x5 2/3")

Weight
2 channels 2.0 kg (5 lb)
4 channels 2.1 kg (5 lb)

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conformément à:
EN 61326-1, EN 61010-1

ISMATEC®

Learn more
about the
Ismatec® product
line by visiting:
www.ismatec.com

For ordering and technical support, please contact:

North America sales@ismatec.com | 1-800-323-4340 | 1-847-549-7600

Europe sales.europe@ismatec.com | +49 (0) 9377 9203-0

 **Cole-Parmer®**