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Operating-
Instructions

Mode d'emploi

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**Digital-pH-Meter
CG 820**

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Digital pH meter with mV-range

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pH metro digital con gama mV

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Temperatur
Temperature
Température
Temperatura

Nullpunkt
Zero point
Zéro
Punto neutro

Abb. 1
Fig. 1
Photo 1
Fig. 1

Funktion
Function
Fonction
Función

Steilheit
Slope
Pente
Gradiente

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OPERATING INSTRUCTIONS

Commissioning

1. The instrument is to be connected to the mains supply.
Mains voltage: 100 ... 240 V ~, 50/60 Hz, no switch over necessary.
2. Connect electrode*.
3. Set switch to position "pH" or "mV".

Calibration of pH meter to the function of electrode

4. Set switch to position "pH".
5. Select two buffer solutions: one with a pH value in the vicinity of the electrode zero point, e.g. pH = 6.87 (according to DIN 19266 and NBS) or pH = 7.00, and the other one with a pH value differing therefrom by about 3 pH units, e.g. pH = 4.01 (according to DIN 19266 and NBS) or pH = 4.00.
6. Equilibrate temperature of buffer solutions and electrode at 25 °C for instance (with a tolerance of ± 0.2 °C if a reproducibility of ± 0.01 pH units is desired). In case a reproducibility of Δ pH = 0.1 units is required, a special temperature equilibration can usually be dispensed with.
7. Set "°C"-knob to buffer solution temperature.
8. Remove cap from refilling orifice of electrode (does not apply to electrodes with gel filling); rinse electrode with distilled water.
9. Immerse electrode in buffer solution with pH value in vicinity of zero point (e.g. pH = 6.87 or pH = 7.00). Set digital display by means of " Δ pH"-knob to value of buffer solution.
10. Rinse electrode with distilled water before immersing electrode in second buffer solution (e.g. pH = 4.01 or pH = 4.00).
11. Set digital display by means of "mV/pH"-knob to pH value of second buffer solution. Unit is thus adapted ("calibrated") to the electrode characteristics. Rinse electrode.

pH measurement

12. Keep switch in "pH"-position.
13. Set "°C" knob to the temperature of the solution to be measured.
14. For measuring the pH value, immerse electrode in solution to be measured. After adequate tempering time, the pH value is read from the display.

* "Electrode" in this context stands for pH combination electrode or for electrode combination, comprising separate indicator electrode and reference electrode.

mV-measurement

15. Commence operation of unit as described under paras 1 to 3.
16. Set switch to "mV"-position.
17. Immerse electrode combination (e.g. platinum-electrode combination for redox potential measurements) in solution to be measured.
18. The measured values are read from the digital display directly.

Note

19. The unit must not be opened on principle for reasons of function.
20. The manufacturer can not accept any warranty claim after unauthorized tampering inside the unit careless use, or deliberate damage to the unit, on the part of the persons handling the unit.
21. Please also observe the instructions for electrodes.

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Warning and safety instructions

For functional and safety reasons, the instrument should be opened only by trained persons. Warranty becomes void in the case of unauthorized tampering inside the instrument or careless or deliberate damage. For safety reasons, the instrument may be used only for the purposes described in this operating instructions manual.

This instrument conforms to protection class I. It has been built and tested in accordance with DIN 57 411, Part 1/VDE 0411, Part 1, Protection measures for electronic measuring equipment and it has left the factory in safe and proper condition. In order to preserve this condition and to ensure safe operation, the user must pay due attention to the instructions and warnings contained in this operating instructions manual.

Before switching on, make sure that the mains voltage setting on the instrument corresponds to the nominal value of the local mains voltage. The voltage setting is specified on the type plate.

Please also observe the corresponding operating instructions for the connected equipment units.

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SCHOTT-GERÄTE GmbH