Leica BM E

Instruction Manual
Gebrauchsanleitung
Mode d'emploi
Manual de instrucciones
Manuale di istruzione

Leica

Leica BM E

Instruction Manual

Contents

1.0 Introduction	1
2.0 Unpacking	1
3.0 Preparation For Use	3.1 Objectives.....	1
	3.1.1 Viewing Bodies.....	1
	3.2 Eyepieces.....	2
	3.3 Substage Condensers	2
4.0 Operation	4.1 Illumination	2
	4.1.1 Built-in Illumination.....	2
	4.1.2 External Illuminator.....	2
	4.2 Stage Preparation.....	2
	4.3 Focusing	2
	4.4 Examining the Specimen at High Power	2
	4.5 Oil Immersion Technique	3
5.0 Care of Microscope	5.1 General	3
	5.2 Mechanical Maintenance	3
	5.3 Lamp Replacement.....	4
	5.4 Electrical Considerations/ Equipment Ratings.....	4
6.0 Fuse Replacement	4
7.0 Anti-Fungus / Mold Warranty	4
8.0 Notes	4

1.0 Introduction

Thank you for purchasing the Leica BM E Compound Microscope. The BM E is designed specifically for student use, with high quality glass optics and a variety of accessories, the BM E Compound Microscope is the ideal choice to meet diverse applications.

2.0 Unpacking

The Leica BM E Compound Microscope is available in preconfigured models. Check the components against the standard equipment list below.

The standard Leica BM E Compound Microscope components are:

1. Stand - includes supporting arm, focusing mechanism, stage, condenser or disc diaphragm, nosepiece and illumination system (mirror or in-base illuminator)
2. Viewing Body - monocular, binocular, or monocular teaching, as ordered
3. Objectives - as ordered
4. Eyepiece(s) - as ordered
5. Immersion Oil - as ordered
6. Dust Cover
7. Allen wrench

Optional accessories such as a carrying case, video and photography attachments are not part of the standard Leica BM E components. If ordered, these items will be shipped separately.

3.0 Preparation for Use

3.1 Objectives

For easy identification, the objectives are color-coded as follows:
4x = red 10x = yellow 40x = dark blue 100x = white

3.1.1 Viewing Bodies

Monocular and Monocular Teaching Bodies: The viewing bodies are secured to the stand by tightening the body thumbscrew on the top right of the instrument then reversing it one-half turn. This permits the viewing body to rotate 360°, but prevents it from being removed.

Binocular Viewing Body: Locate the knurled screw on the top right of the microscope stand and loosen it. Rotate the viewing body to the desired position. Tighten the knurled screw firmly but without overtightening. The viewing body is now precisely located on the optical axis regardless of the selected orientation.

Explanation of Warning Symbols

CAUTION: Replace with same type and rating fuse (IEC 60127-2/SS 3 standard).

CAUTION: Risk of electric shock

CAUTION: (Refer to accompanying documents)

3.2 Eyepiece

The eyepiece is secured to the eyepiece tube. Binocular viewing bodies do not utilize set screws to secure the eyepieces.

3.3 Substage Condenser

The Leica BM E Compound Microscope is equipped with one of two available substage condensers:

- 0.65 N.A. condenser, a five-position aperture disc diaphragm which is color coded to match the objective markings.
- 1.25 N.A. (fixed) Abbe Condenser on a spiral mount

4.0 Operation

CAUTION: *The top collector lens of the illumination system in the base of the microscope may become warm or hot. Avoid touching the lens directly.*

The Leica BM E microscope should always be used on a hard, stable and level surface such as a laboratory table or work bench.

4.1 Illumination

4.1.1 Built-in Illumination System

The Leica BM E Compound Microscope should only be plugged into a suitably grounded electrical outlet.

To turn on the electrical power, flip the switch located on the back left of the base of the BM E. If you have the Tungsten-Halogen model, the intensity control which is located on the left side of the base should be turned all the way down before turning the instrument off.

4.1.2 External Illuminator and Mirror

Position the external illuminator approximately 5 cm (2") in front of the mirror.

NOTE: *Not all Leica BM E Compound Microscopes are equipped with a mirror. While looking through the eyepiece adjust the mirror to provide the brightest and most evenly illuminated field-of-view.*

4.2 Stage Preparation

Lower power objectives have a greater depth of field and are generally used for initial focusing and viewing. The lower the objective power, the greater area of specimen surface will be seen in the field of view.

Before placing a specimen on the stage, rotate the nosepiece until the lowest power objective is in the viewing position.

Be certain that the stage surface is free of dust, grit, or any other material. These will scratch the slide or stage or may interfere with the movement of a slide across the surface of the stage.

Place the slide under the stage clips with the cover glass facing upward. If using a graduated mechanical stage, place the slide between the stage fingers.

Position the specimen area of the slide over the center of the stage aperture. When using a graduated mechanical stage, use the stage control knobs to move the specimen slide.

Raise the stage by turning the coarse focus knob until the stage reaches the upward stop. Do not view the slide through the eyepiece at this point.

The Leica BM E Compound Microscope is ready to be focused.

4.3 Focusing

Using the 4x or 10x objective, focus the Leica BM E Compound Microscope by looking through the eyepiece and turning the fine focus knob. Continue to view the specimen through the eyepiece until a crisp, clear image appears.

Adjust the aperture disc diaphragm condenser to match the objective being used to obtain the clearest possible image. Adjust the fine focus knob to sharpen the image to the center of the field of view.

If using the 1.25 N.A. Abbe condenser, focus by moving it up and down in its spiral sleeve to the point at which the field is most evenly illuminated. Adjust the iris diaphragm by sliding the control lever to obtain the clearest image.

The position of the aperture must be reset when changing to higher power objectives. As magnification increases, the aperture diaphragm must be opened.

4.4 Examining Specimen at High Power

When examining a specimen at high magnification, position the slide so the feature is centered in the field of view. If the feature you want to observe appears to the left of the field of view, move the slide left to center it. If the feature you want to observe appears to the right of the field of view, move the slide to the right.

By rotating the nosepiece, bring the next higher power objective into viewing position. The objective is properly aligned when the nosepiece clicks into place.

4.5 Oil Immersion Technique

To utilize the full numerical aperture of an immersion objective (with inscription "oil"), oil is added to both the objective and specimen using the following procedure:

1. Focus on the specimen with a lower power objective (40x).
 2. Put a drop of immersion oil on the specimen slide.
 3. Turn the nosepiece to bring the immersion objective into the light path, and focus with the fine adjustment knobs.
- Care should be taken to prevent bubbles from forming in the oil, since they will deteriorate the lens performance. Simply reapply the immersion oil if bubbles form.

IMPORTANT: *After each use of the 100x objective, use a lens tissue or soft cloth to wipe off all traces of oil from the objective and the slide.*

5.0 Care of the Microscope

5.1 General

It is important to keep all optical components clean to obtain the best optical performance. Keep a dust cover on the microscope when it is not in use.

The objectives are easily subjected to dirt, dust and oil. Whenever lack of contrast, cloudiness or poor definition is encountered, carefully check the condition of the front objective lens with a magnifier. To assure all objective lenses are clean, periodically examine them with a magnifying glass.

If an optical surface becomes badly coated with dust or dirt, use an air-filled syringe or soft brush to remove it before wiping the surface clean.

IMPORTANT: *Clean optical surfaces with a lint-free cloth, lens tissue or a cotton swab moistened with methanol or commercially available glass cleaner. Avoid excessive use of solvents. Interior surfaces can be loosened or badly damaged if solvents seep in around the lenses.*

The higher magnification objectives use a small front lens. The surface of this front lens can be cleaned with a small cotton swab. Dampen the cotton tip with methanol and squeeze it dry. Wipe the front lens lightly. Make sure that the cotton tip contacts the lens surface. Check the objective with a magnifier after the cleaning.

IMPORTANT: *If it is necessary to remove the viewing body, DO NOT TOUCH the outer surface on the underside of the viewing body. Finger prints on this surface will reduce the image clarity. This lens can be cleaned in the same manner as objectives and eyepieces.*

5.2 Mechanical Maintenance

All non-optical components of the Leica BM E Compound Microscope can be cleaned with soap and water. Avoid using all other harsh solvents or acetone-based products when cleaning.

The Leica BM E Compound Microscope requires only routine maintenance. Neither the focusing mechanism nor the mechanical stage require lubrication. Contact your Leica customer service representative for more information.

5.3 Lamp Replacement

 The Leica BM E utilizes a 7 watt fluorescent lamp with a life expectancy of 10,000 hours or a 12 watt tungsten-halogen lamp with a life expectancy of 2,000 hours. When replacement is necessary contact a qualified service technician.

 **CAUTION:** The glass envelope of the lamp may be extremely hot. Do not attempt to change the lamp before it is completely cooled or without wearing adequate skin protection. The glass should not be handled with bare hands. If the lamp is touched with bare hands clean the surface with a lens tissue and methanol prior to installation.

5.4 Electrical Considerations / Equipment Ratings

 **CAUTION:** As a safety precaution, the power cord has been grounded to the microscope base. Never use an adapter between the power cord and the power source; it will render the grounding feature ineffective.

 **CAUTION:** Risk of electric shock. Removal of the bottom cover of the microscope exposes hazardous voltages, which if contacted, could cause injury or death. Please refer servicing to a qualified service technician.

Supply Voltage and Range Supply Frequency

Tungsten Halogen models:

- 120VAC $\pm 10\%$ - 50/60Hz, input power - 15VA
- 220-240VAC $\pm 10\%$ - 50/60Hz, input power - 15VA

Fluorescent models:


- 120VAC $\pm 10\%$ - 60Hz, input power - 20VA
- 220-240VAC $\pm 10\%$ - 50/60Hz, input power - 45VA

Range of Environmental Conditions

- Indoor use only
- Temperatures 5°C to 40°C (41°F to 104°F)
- Maximum Relative Humidity - 80% for temperatures up to 31°C (87.8°F) decreasing linearly to 50% relative humidity at 40°C (104°F)
- Installation Category II (Overvoltage Category)
- Pollution degree 2
- Altitude up to 2,000 meters

6.0 Fuse Replacement

 **CAUTION:** The fuse can be replaced only with the same size and rating (IEC 60127-2/SS 3 standard).

Tungsten Halogen models: 

- 120 VAC - 200 mA/250VAC, Time-Lag
- 220-240 VAC - 160 mA/250VAC, Time-Lag

Fluorescent models: 

- 120 VAC - 250 mA/250VAC, Time-Lag
- 220-240 VAC - 250 mA/250VAC, Time-Lag

7.0 Anti-Fungus / Mold Warranty

Fungus and mold are more likely to grow on optical surfaces in high temperature and high humidity regions. The standard anti-reflective optical coatings used on the DM E, CM E and BM E microscopes protect the optical surfaces from these climatic conditions. Microscopy products from Leica Microsystems have a history of performing well for long periods of time without fungus growth. This confidence allows us to include growth of fungus in your standard Warranty for DM E, CM E and BM E.

8.0 Notes

Leica assures the quality of our instruments to be free of all faults in material and manufacturing. We do not, however, cover damage resulting from improper or careless handling. Contact your Leica representative for more information.

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Leica Microsystems Inc.