Leica CM E

Instruction Manual
Gebrauchsanleitung
Mode d'emploi
Manual de instrucciones
Manuale di istruzione

Leica

Leica CM E

Instruction Manual

Contents

1.0 Introduction2
2.0 Unpacking2
3.0 Preparation For Use	
3.1 Viewing Body.....	2
3.2 Installing Objectives.....	3
3.3 Substage Condenser	3
3.4 Filters	3
4.0 Operation3
5.0 Care of Microscope	
5.1 General.....	3
5.2 Mechanical Maintenance.....	4
5.3 Lamp Replacement	4
5.4 Electrical Considerations/ Equipment Ratings	4
6.0 Fuse Replacement5
7.0 Anti-Fungus / Mold Warranty5
8.0 Notes5

1.0 Introduction

Thank you for purchasing the CM E from Leica Microsystems. The Leica CM E Compound Microscope combines exclusive design features and a full range of accessories making it a truly versatile, quality instrument for advanced academic and routine clinical applications.

CAUTION: This instrument should only be used as described. Hazards may result if used improperly.

2.0 Unpacking

Carefully remove the microscope components from the packing carton. Check the components against the following standard equipment list:

1. Stand: includes focusing mechanism, graduated mechanical stage, nosepiece and integral illumination system.
2. Viewing body
3. Objectives
4. Eyepiece(s)
5. Condenser
6. Immersion oil
8. Dust cover
9. Power cord (as ordered)
10. Blue filter
11. 3mm Allen wrench

Optional items such as contrast accessories, camera adapters and systems and carrying cases are not shipped as part of the standard equipment. These items are delivered in separate packages.


Please do not discard any of the packing materials. They should be used for safely storing and transporting the instrument should the need arise.

3.0 Preparation for Use

3.1 Viewing Body

Attach the viewing body to the stand by backing off the set screw on top of the stand with the tool provided. Then insert the dovetail of the viewing body into the stand support and tighten the set screw firmly without forcing it. This precisely locates the body on the optical axis of the microscope, regardless of the selected orientation. After the body is attached, insert the eyepieces into the eyetubes. Attach the power cord and turn on the CM E using the switch on the back of the instrument. The knurled dial on the lower left side will control the illumination intensity.

Explanation of Warning Symbols

	CAUTION: Replace with same type and rating fuse (IEC 60127-2/SS 3 standard).
	CAUTION: Risk of electric shock
	CAUTION: (Refer to accompanying documents)

Interpupillary distance is changed by simply spreading or pushing eyetubes on the binocular and trinocular bodies.

If you have purchased a standard or custom outfit, you will notice the objectives are already installed on the nosepiece and the condenser is in place under the stage. In this case, go to the OPERATION section. If you purchased your CM E in components, please continue.

3.2 Installing Objectives

Attach the objectives to the revolving nosepiece holes counterclockwise from the lowest magnification to the highest. *Always rotate the nosepiece by using the knurled nosepiece ring.*

3.3 Substage Condenser

Adjust the condenser mount to its lowest position, then insert the substage condenser into the condenser mount underneath the stage. Tighten the condenser thumbscrew. Each condenser is furnished with an iris diaphragm which controls the effective numerical aperture of the objective. To open and close this diaphragm, simply slide the condenser control to the right or to the left.

3.4 Filters

The blue filter is factory installed in the underside of the condenser. For phase contrast applications replace the blue filter with the green filter. Special care must be used in handling filters. Scratches and fingerprints will result in poor images and lower light levels.

4.0 Operation

CAUTION: As a safety precaution, the power cord has been grounded to the microscope base. Never use an adapter between the power cord and the power source; it will render the grounding feature ineffective.

CAUTION : The top collector lens of the illumination system in the base of the microscope may become warm or hot. Avoid touching the lens directly.

4.1 Always use your microscope on a hard, stable surface.

4.2 Plug the microscope power cord into a suitable grounded electrical outlet. A grounded 3-wire cord is provided.

4.3 Turn on the microscope with the switch in back and rotate the illumination control knob on the bottom left side of the instrument.

4.4 Set the illumination control knob to the lowest setting. The illumination control allows you to adjust the intensity of light produced by the illumination system.

4.5 Fully open the aperture diaphragm of the condenser by sliding the control to the extreme left.

4.6 Using the substage condenser focusing knob, raise the condenser to the top of its excursion.

4.7 Place a specimen slide on the stage.

4.8 Rotate the nosepiece to move the 4x objective into working position.

4.9 Raise the stage by rotating the coarse adjustment knob to its positive stop. Using the fine adjustment knob, bring the specimen into sharp focus. The CM E has been factory adjusted such that stage focus can be found ≤ 1.5 rotations of highest focus adjustment.

4.10 Adjust the eyetubes for interpupillary distance (IPD) and eye difference, then set the IPD # for both eyepiece tubes to match the IPD number on the binocular viewing body scale located between the eyepiece tubes. Both eyepiece tubes are focusable to compensate for refractive differences of the eyes. To correctly set the eyetubes, rotate the 40x objective into place and focus on the specimen. Then rotate the 4x objective into place and focus each eyetube separately using the focus ring on each eyetube. Do not use the coarse or fine focus controls on the side of the instrument at this point.

4.11 Remove an eyepiece and view the back aperture of the objective. Close the condenser aperture diaphragm and then, to obtain the full resolving power of the microscope, reopen until the diaphragm leaves just disappear from view. Replace the eyepiece. The aperture diaphragm can be adjusted to enhance contrast and/or increase the depth of focus.

4.12 When changing objectives, the positions of the aperture diaphragm must be reset. As magnification increases, the aperture diaphragm must be opened as required.

5.0 Care of the Microscope

5.1 General

Keeping all optical components clean is important for maintaining good optical performance. The microscope should always be covered with the plastic dust cover provided with the instrument when it is not in use. If any optical surface becomes badly coated with dust or dirt, it should be cleaned by blowing it off with a syringe or brushing it off

with a camel hair brush before attempting to wipe the surface clean.

Optical surfaces should be cleaned with a lint-free cloth, lens tissue or a cotton swab moistened with methanol or commercially available glass cleaner. It is very important to avoid the excessive use of solvent.

The cloth lens tissue or cotton swab should be moistened with solvent, but not be wet enough for the solvent to seep around the lens.

No part of the microscope is quite so vulnerable to dirt, dust and oil as the front lens of the objective. Whenever lack of contrast, cloudiness or poor definition is encountered, carefully check the condition of the front lens with a magnifier.

The lower magnification objectives have fairly large front lenses and can be cleaned with a cloth or lens tissue wrapped around a finger and moistened with methanol. The 40x and 100x objectives require a little more care and examination with a magnifier. To achieve the high degree of flatness obtained with the higher magnification objectives, it is necessary to use a small concave front lens of fairly short radius or curvature. The surface of this front lens can be readily cleaned with a toothpick covered with a cotton tip, or with a small cotton swab. Moisten the cotton with methanol and squeeze it almost dry. Wipe the front lens lightly without applying undue force or scrubbing action. Make sure that the cotton tip contacts the concave lens surface. Check the objective with a magnifier after cleaning.

Should you need to remove the viewing body, be careful not to accidentally touch the outer lens surface on the underside of the body. Finger prints on this surface will reduce image clarity. This lens can be cleaned in the same manner as objectives and eyepieces.

5.2 Mechanical Maintenance

Using the dust cover is the most important factor in keeping the microscope in good mechanical and physical condition.

The Leica CM E has a durable finish. All non-optical parts of the instrument can be cleaned with methanol, naphtha or soap and water. Avoid the use of all other organic solvents when cleaning the instrument.

5.3 Lamp Replacement

CAUTION: Unplug the instrument before attempting to change the lamp.

CAUTION: The glass envelope of the lamp may be extremely hot. Do not attempt to change the lamp before it is completely cooled or without wearing adequate skin protection. The glass should not be handled with bare hands. If the lamp is touched with bare hands clean the surface with a lens tissue and methanol prior to installation.

The only component in the Leica CM E that may require periodic replacement is the illuminator lamp.

12V, 20W Halogen replacement lamp # 13496025

- A. Before replacement, unplug the instrument.
- B. Open the lamp door by pulling the locking pin
- C. After it has cooled, carefully hold the lamp using a lens cloth and remove it from its socket.
- D. Place the new lamp into the socket using a lens cloth.
- E. Close the lamp door, being sure to push the locking pin in until it stops.
- F. Never operate the microscope illuminator unless the lamp door is securely closed.

5.4 Electrical Considerations / Equipment Ratings

CAUTION: Risk of electric shock. Removal of the bottom cover of the microscope exposes hazardous voltages, which if contacted, could cause injury or death. Please refer servicing to qualified personnel.

Supply Voltage and Range Supply Frequency

Tungsten Halogen models:

- 100 VAC $\pm 10\%$ - 50/60 Hz, input power - 28 VA
- 120 VAC $\pm 10\%$ - 50/60 Hz, input power - 28 VA
- 220-240 VAC $\pm 10\%$ - 50/60 Hz, input power - 28 VA

Range of Environmental Conditions

- Altitude up to 2000m
- Indoor use only
- Temperatures 5°C to 40°C
- Installation Category II (Overvoltage Category)
- Pollution degree 2

6.0 Fuse Replacement

CAUTION: Unplug the instrument before changing any fuses.

CAUTION: The fuse can be replaced only with the same size and rating (IEC 60127-2/SS 3 standard).

- 100 VAC - 400mA/250VAC, Time-Lag
- 120 VAC - 315mA/250VAC, Time-Lag
- 220-240 VAC - 200mA/250VAC, Time-Lag

7.0 Anti-Fungus / Mold Warranty

Fungus and mold are more likely to grow on optical surfaces in high temperature and high humidity regions. The standard anti-reflective optical coatings used on the DM E, CM E and BM E microscopes protect the optical surfaces from these climatic conditions. Microscopy products from Leica Microsystems have a history of performing well for long periods of time without fungus growth. This confidence allows us to include growth of fungus in your standard Warranty for DM E, CM E and BM E.

8.0 Notes

Leica assures the quality of our instruments to be free of all faults in material and manufacturing. We do not, however, cover damage resulting from improper or careless handling. Contact your Leica representative for more information.

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 5x20 mm, 3,15 A/250V, ad azione rapida (n° 13WFAG30001)

Due to a policy of continuous development, we reserve the right to change specifications without notice.

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Tel 716 686 3000
Fax 716 686 3085
www.leica-microsystems.com
ISO-9001 Certified

The Leica logo is written in a classic, elegant script font. The letters are black and have a slight shadow or depth, giving it a three-dimensional appearance. The 'L' is particularly large and stylized, with a long horizontal stroke that curves under the 'e'. The 'i' has a small dot, and the 'c' is a simple, rounded shape. The overall look is sophisticated and professional.

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Leica Microsystems Inc.