




**Micromanipulator System**  
**NT-88NE**

**Water-Hydraulic  
Micromanipulators**  
**MW0-3/MHW-3**



## Developed Using the Latest Micromanipulation Technology, the NT-88NE and Related Equipment Combine Outstanding Precision with Easy Operation

The NT-88NE standard package consists of:

3D Motor-Driven Coarse Manipulator MM-188NE: 2 pcs

3D Hydraulic Fine Micromanipulator MO-188NE: 2 pcs

Universal Joint U-NE: 2 pcs

Pipette Holder HI-7 with Rotating Mechanism HIR: 2pcs

Holder Keeper HIK-3: 2 pcs

Adapter for TE2000 NS-888NE: 1pc

To this equipment, add your desired injectors to build your system.

Whatever the application—IVF (in-vitro fertilization), ICSI (intracytoplasmic sperm injection), electrophysiology, or biotechnology—the NT-88NE and other top-notch micromanipulation equipment from Nikon will meet the task.



## 3D Hydraulic Fine Micromanipulator MO-188NE



Hanging-type joystick for MO-188NE

### Universal Joint "U-NE"

This convenient mechanism simplifies operation by allowing the pipette holder to be moved in a single direction. Angle adjustment for both vertical and horizontal directions is possible via two separate control knobs. This product can be purchased individually.



Drive unit	
X-, Y-, Z-axes movement	10 mm (max.)
Control knob	
X-, Y-, Z-axes control knob movement	10 mm (max.)
Joystick movement (X- and Y-axes)	2 mm (max.)
Control knob movement (one rotation)	250 μm (with 2.5 μm min. scale)
Configurations	Drive unit, control unit (with magnetic base), pipette holder, universal joint, metal plate for magnetic base

X-axis control knobs are located symmetrically on both sides of the joystick to eliminate left side/right side positioning problems. This design allows faster assembly and saves on bench space. Furthermore, a 10mm X-Y operating range provides more versatility at lower magnifications.

## 3D Motor-Driven Coarse Manipulator MM-188NE



Joystick control for MM-188NE

A "Resume" feature returns the pipette to the original position after the pipette has been raised for quicker, easier petri dish changes. The drive unit has individual sets of index calibrations for the X, Y, Z axes for fast repositioning. In addition, a microscope adapter with enhanced resistance to vibrations ensures more stable manipulation, while the joystick's compact, space-saving design and height adjustment provide comfortable operation.

Drive unit	
X-, Y-, Z-axes movement	22 mm (max.)
Control unit	
Minimum speed	Approx. 50 μm/sec.
Maximum speed	Approx. 250 μm/sec.
High-speed mode speed	1.6 mm/sec.
Joystick movement	X-, Y-axes
Up-down button movement	Z-axis
High/low speed adjustment	Provided
Configurations	Drive unit, control unit, power transformer

## 3D Manual Coarse Manipulator MN-188NE



A manual 3D coarse positioning manipulator allows you to quickly move the micro-pipette to the ideal position. The MN-188-NE's compact size saves installation space, while a stopper on the Z-axis further improves ease of operation.

Drive unit	
X-, Y-, Z-axes movement	30 mm (X and Y), 30 mm (Z) with stopper
Configurations	Main unit

## Oil-type Microinjector IM-9B / Pneumatic Injector IM-9C



IM-9B

IM-9C

These injectors provide smoother and more stable maneuverability by applying a new mechanism to minimize backlashes. Also, its 53mm working distance, large-volume syringe and just-the-right-size control knob with graduated scales all contribute to easy control of the injection. The IM-9B has a finely-ground metal syringe for smoother movement, while its new tube connector with a multi-purpose valve makes for easier filling of oil without generating air bubbles.

	IM-9B	IM-9C
Movement distance		
Movement range	53mm	53mm
Full rotation of knob	500μm	6mm
Minimum graduation	10μm	120μm
Control value		
Full rotation of knob	10μl	480μl
Minimum graduation	0.2μl	9.6μl
Syringe	1.1cc metal type	4.2cc glass type

### Pipette Holder HI-7 with Rotating Mechanism HIR



This is a standard accessory included with these microinjectors. Featuring a rotating mechanism, the HIR, this accessory enhances the easy alignment of the pipette tip in the microscope field, while reducing oil leaks.

## 1D Hydraulic Micromanipulator MO-109NE



The 1D micromanipulator is attached to a ball-joint, enabling versatile angle adjustment.

It can also be used in combination with other systems for very effective injection. Attaching this equipment to the MO-188NE creates the 4D Hydraulic Micro-manipulator MO-189NE.

4D Hydraulic Micromanipulator MO-189NE



Drive unit	
Movement in one direction (T-unit)	10 mm (max.)
Control unit (T unit)	
T-unit movement	10 mm (max.)
Control knob movement (one rotation)	250 μm (with 2.5 μm min. scale)
Configurations	Drive unit (with T-unit ball joint), control unit (with mechanical T-handle)



## Two-Electrode Holder HD-21

The HD-21 allows you to install two pipette holders close to each other on a 3D hydraulic micromanipulator. It's ideally suited for embryo-biopsy applications with controls to adjust angles smoothly, thanks to the ball-joint pivot construction.



## Combination Examples

**Eclipse TE2000 with  
MM-188NE x 2,  
MO-188NE x 2,  
IM-9B, IM-9C**




**TS100 with MM-188NE x 2,  
MO-188NE x 2,  
IM-88, IM-88H**



**Stereoscopic Microscope  
SMZ1500 with  
MM-188NE, MO-188NE,  
IM-9B, IM-9C**

## System Chart




## NT-88NE Combination Chart

	Embryo-biopsy set	In-vitro fertilization set	Biotechnology set	Electrophysiology set*
4D Hydraulic Fine Micromanipulator MO-189NE with hanging-type joystick control	-	-	1 pc. (For injection side)	-
3D Hydraulic Fine Micromanipulator MO-188NE with hanging-type joystick control	2 pcs.	2 pcs	1 pc. (For holding side)	-
3D Water-Hydraulic Micromanipulators MHW-3/MWO-3	-	-	-	2 pcs. of MHW-3 or MWO-3
3D Motor-Driven Coarse Manipulator MM-188NE	2 pcs.	2 pcs.	2 pcs.	-
3D Manual Coarse Manipulator MN-188NE	2 pcs. (For economical set)	2 pcs. (For economical set)	2 pcs. (For economical set)	-
1D Hydraulic Fine Micromanipulator MO-109NE	-	-	-	-
Pneumatic Microinjector IM-9C	1 pc. (For holding side; economical set)	1 pc. (For holding side; economical set)	1 pc. (For holding side; economical set)	-
Oil-type Microinjector IM-9B	2 pcs.	2 pcs.	2 pcs.	-
Two-Electrode Holder HD-21	1 pc.			
Ball-Joint B-8C/Universal Joint "U-NE"	-	-	B-8C 2 pcs.	-
Pipette Holder HI-7	2 pcs.	2 pcs.	2 pcs.	-

\*Other items for electrophysiology are available from other manufacturers. Mounting adapters are available to fit to your microscope.

## Optional Equipment for Making Customized Pipettes

### Pipette Grinder EG-400



The EG-400 virtually eliminates vibration and drift over the grinding surface due to a specially designed micro-motor on the variable-speed rotating wheel.

Grinding wheel	Diamond
Adjustment	Coaxial coarse/fine manipulator
Accessories	Pipette holder, motor output adjustment dial, cleaner

### Microforge MF-900



The MF-900 is designed to produce custom-made injection and holding pipettes.

Microscope	Eyepieces (10X), objectives (5X, 10X), trans-illumination
Accessories	Heater, heater output adjustment dial, heater foot switch, air hose

### Glass Microelectrode Puller PN-30



The PN-30 is designed to produce microelectrodes from 1 mmø to 1.5 mmø glass capillaries.

Heater voltage	2V
Current	45A
Glass capillaries	1 - 1.5 mmø

	Embryo-biopsy set	In-vitro fertilization set	Biotechnology set	Electro-physiology set
Micropipette GD-1/G-100	G-100 1 pc.	G-100 1 pc.	GD-1/G-100 1 pc. each	GD-1 1 pc.
Micropipette Puller PC-10/PN-30	PN-30 1 unit	PN-30 1 unit	PC-10 1 unit	PC-10 1 unit
Pipette Grinder EG-400/EG-6	EG-400 1 unit	EG-400 1 unit	EG-400 1 unit	EG-6 1 unit
Microforge MF-900	1 unit	1 unit	1 unit	1 unit

## 3D Water-Hydraulic Micromanipulator MWO-3

The MWO-3 features a cutting-edge, hydraulic double-cartridge mechanism to minimize drifts caused by temperature changes. This design ensures injections with a constant pressure for stable operation that is free from drift, vibration, and electromagnetic disturbances.

Movement range	Fine 2mm (X, Y and Z axes), coarse 30mm (X, Y and Z axes)
Full rotation of knob	Fine 50 $\mu$ m, coarse approx. 4mm
Minimum graduation	Fine 0.2 $\mu$ m



MWO-3

## 3D Water-Hydraulic Micromanipulator MHW-3

The MHW-3 features minimized drift due to temperature changes, a scale increment as small as 0.2 $\mu$ m, and a fine control manipulation with a one-to-five moving ratio. Also, because Teflon tubes are secured to the microscope main body, vibration of the tip of the micropipettes has been eliminated, resulting in stable micromanipulation.

Movement range	Fine 2mm (X, Y and Z axes), coarse 30mm (X, Y and Z axes)
Full rotation of knob	50 $\mu$ m
Minimum graduation	Fine 0.2 $\mu$ m



MHW-3 with NN-L adapter

## Manual 3D Micromanipulators NMN-21/25

The NMN-21 and -25 are designed for drift-free operation, making them ideal for delicate, patch-clamp research. The unit has no hydraulic system nor spring mechanism and its unique double-slider structure minimizes vibrations transmitted from the operator's hand to the tip of the pipette via the control knob. The unit's motor-less and actuator-less design also eliminates electromagnetic disturbances. The NMN-25 features a compact, space-saving design without a scale.

	NMN-21	NMN-25
Movement range	Fine 6mm (X, Y and Z axes), coarse 15mm (X, Y and Z axes)	Fine 6mm (X axis), 10mm (Y axis), 6mm (Z axis); coarse 15mm (X axis), 10mm (Y axis), 15mm (Z axis)
Full rotation of knob	Fine 250 $\mu$ m, coarse approx. 4mm	Fine 250 $\mu$ m (X and Z axes), 292 $\mu$ m (Y axis)
Minimum graduation	Fine 1 $\mu$ m	—



NMN-21

**WARNING** TO ENSURE CORRECT USAGE, READ THE CORRESPONDING MANUALS CAREFULLY BEFORE USING YOUR EQUIPMENT.

Specifications and equipment are subject to change without any notice or obligation on the part of the manufacturer. December 2004  
©1997-2004 NIKON CORPORATION

Micromanipulators and related equipment other than microscopes are manufactured by:

**NARISHIGE**  
NARISHIGE SCIENTIFIC INSTRUMENT LAB.  
27-9, Minamikarasuyama 4-chome, Setagaya-ku, Tokyo 157-0062, Japan  
Phone: +81-3-3308-8233 Fax: +81-3-3308-2005



ISO 14001 Certified  
NIKON INSTRTECH CO., LTD.



ISO 9001 Certified  
NIKON CORPORATION  
Instruments Company



ISO 14001 Certified  
NIKON CORPORATION  
Yokohama Plant

### NIKON INSTRTECH CO., LTD.

Parale Mitsui Bldg., 8, Higashida-cho, Kawasaki-ku,  
Kawasaki, Kanagawa 210-0005, Japan  
phone: +81-44-223-2175 fax: +81-44-223-2182  
<http://www.nikon-instruments.jp/eng/>

### NIKON INSTRUMENTS (SHANGHAI) CO., LTD.

CHINA phone: +86-21-5058-5055 fax: +86-21-5058-5060  
(Beijing office)

CHINA phone: +86-10-5869-2255 fax: +86-10-5869-2277

### NIKON SINGAPORE PTE LTD

SINGAPORE phone: +65-6559-3618 fax: +65-6559-3668

### NIKON MALAYSIA SDN. BHD.

MALAYSIA phone: +60-3-78763887 fax: +60-3-78763387

### NIKON INSTRUMENTS EUROPE B.V.

P.O. Box 222, 1170 AE Badhoevedorp, The Netherlands  
phone: +31-20-44-96-222 fax: +31-20-44-96-298  
<http://www.nikon-instruments.com/>

### NIKON FRANCE S.A.S.

FRANCE phone: +33-1-45-16-45-16 fax: +33-1-45-16-00-33

### NIKON GMBH

GERMANY phone: +49-211-9414-0 fax: +49-211-9414-322

### NIKON INSTRUMENTS S.p.A.

ITALY phone: +39-55-3009601 fax: +39-55-300993

### NIKON AG

SWITZERLAND phone: +41-43-277-2860 fax: +41-43-277-2861

### NIKON UK LTD.

UNITED KINGDOM phone: +44-20-8541-4440 fax: +44-20-8541-4584

### NIKON INSTRUMENTS INC.

1300 Walt Whitman Road, Melville, N.Y. 11747-3064, U.S.A.  
phone: +1-631-547-8500; +1-800-52-NIKON (within the U.S.A. only) fax: +1-631-547-0306  
<http://www.nikonusa.com/>

### NIKON CANADA INC.

CANADA phone: +1-905-625-9910 fax: +1-905-625-0103



### NIKON CORPORATION

Fuji Bldg., 2-3, Marunouchi 3-chome, Chiyoda-ku,  
Tokyo 100-8331, Japan  
[www.nikon.com/](http://www.nikon.com/)