## Biotest HYCON

RCS High Flow  
Luftkeimsammler/Air Sampler

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
### Table of Contents

		Page
1.	Individual Parts and Functions	3, 49
2.	Application	26
3.	Principles of Operation and Construction	27
4.	Operation Instructions	28
4.1	General Notes on Operation	28
4.2	Sterilization and Decontamination Before Use	29
4.3	Insertion of the Agar Strip	29
4.4	Switching the Instrument On	30
4.5	Setting the Sample Volume	30
4.6	Starting the Instrument	31
4.7	Switching the Instrument Off	31
4.8	Removal of the Agar Strip	31
4.9	Setting the Delay Time	32
4.10	Power Supply	32
4.11	Tripod and Table-Top Tripod	32
4.12	RS232-Interface	33
5.	Evaluation of Results	33
6.	Service-Menu	34
6.1	Bluetooth- Data transfer	34
6.2	HYCON-ID	34
6.3	Delay time	35
6.4	Individual sampling volume	35
6.5	Interval-Sampling	36
6.6	Setting Time	37
7.	Display Messages/ Acoustic Signals	38
8.	Service and Maintenance	40
8.1	General RESET	40
8.2	Calibration	40
8.2.1	Electronic calibration	42
8.3	Cleaning of the Rotor	42
9.	Accessories	43
10.	Technical Data	45
11.	Warranty	46
12.	Waste Disposal	46

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## 2. Application

The Biotest RCS High Flow Air Sampler is used to investigate the microbiological quality of ambient air, the functionality of air treatment equipment, and the effectiveness of decontamination measures. It is used routinely in areas where specified airborne microbial counts, as described by law or through individual standards, may not be exceeded.

The RCS High Flow Air Sampler provides a simple means of testing concentrations of airborne microbes in the following areas:

- Sterile and aseptic production lines in pharmaceutical industry (Cleanroom and Isolator environments)
- Surgery theatres and intensive care units in hospitals
- Food and beverage industry
- Cosmetic industry
- Packaging industry
- General indoor and outdoor air quality sampling

## 3. Principles of Operation and Construction

The Biotest RCS High Flow Air Sampler operates on the impaction principle. It facilitates the quantitative separation of airborne microorganisms from sample volumes ranging from 10 to 1500 liters.

The air stream enters the rotor at the front of the instrument, is set into rotation by the movement of the fan, and separates microbes contained in the air onto the agar strip through centrifugal force. To prevent the occurrence of turbulence in the air inlet channel, the air outlet passage is directed to the rear, parallel with the instrument. The air flows at a constant rate which is approximately 100l/min at a rotor speed of 8200 rpm. The instrument can be calibrated to ensure that the sampling volume remains exact.

The instrument is constructed with special attention to ergonomics. It operates independently of a mains power supply using a 9.6 Volt rechargeable nickel cadmium battery. The battery is easy to change through a door on the battery compartment at the rear of the unit. The RCS High Flow can also be operated with an optional mains power supply.

Built-in electronic circuitry controls the speed and the running time of the motor and simultaneously monitors the battery voltage. Sampling volumes of 10, 20, 50, 100, 200, 500 and 1000 liters have been permanently programmed into the instrument. In addition, it is possible to enter three individually chosen sampling volumes of up to a maximum of 1999 liters

(it is recommended not to exceed 1000 l) as well as a Delay Time in the range of 1 minute to 1 hour and 59 minutes. Operation of the instrument is via a keyboard with integrated display which is in front of the carrying handle, allowing it to be operated with one hand. The instrument provides various memory functions, giving the user a wide variety of important control and security functions.

The rotor, which contains the agar strip, is driven via a magnetic coupling flange. This enables it to be simply removed from the instrument. For safety reasons, a stainless steel protection cap is situated over the rotor, without which the instrument cannot function.

The instrument can be operated in a vertical or a horizontal position but should be protected against strong movements. For sampling in areas difficult to access, the instrument can be attached to a tripod by means of a threaded attachment in the base of the housing.

## 4. Operation Instructions

### 4.1. General Notes on Operation

**Note:** The Biotest RCS High Flow Air Sampler can only be operated with rotors whose serial number starts with the letter »H«.

The standard equipment of the Biotest RCS High Flow Air Sampler includes two rechargeable batteries, one of which is already installed in the instrument, a serial cable for connection with a PC, CD-Rom with RCS Management Software, manual and calibration certificate. For transport it is recommended that the rotor be removed from the instrument and placed in the appropriate compartment in the case. Before using the Biotest RCS High Flow Air Sampler for the first time, the batteries must be charged (see 4.10.). Sterile gloves should be worn when handling the air sampler and agar strips to prevent contamination.

#### Set up of RCS High Flow prior to use

The RCS High Flow is programmed to be in »stand-by« mode on delivery, which means that the expiration date for the calibration of the delivered rotors is not activated.

**Please note:** The first switching on of the RCS High Flow activates the internal clock which in turn calculates the expiration date for the rotor calibration.

A period of 365 days is pre-selected, starting with the day of first use of the instrument. This can be altered with an electronic calibration.

#### Additional Functions of the RCS High Flow Rotor Selection

The instrument is able to store data for 10 different rotors including the respective calibration factors. Each rotor number also includes the calibration expiration date which is shown on the display. After switching the instrument on, the display of available rotor numbers can be selected. By then pressing the [+] and [-] key all available rotor serial numbers can be displayed. By pressing [MENU/ON] a rotor can be selected. The rotor serial number is displayed with the calibration expiration date. If the expiration date is exceeded, the instrument warns with an error message. The alert must be confirmed by pressing the [MENU/ON] key. The alert is then extinguished and in starting mode (QA-Level 1 of RCS Management Software) volume selection is allowed. If a higher QA-Level is selected, the display will show the rotor selection menu again. In this case the volume selection can be performed if a rotor with a valid calibration expiration date was selected.

#### Time-Setting:

Time can be changed through the Service-Menu.

#### Interval-Sampling:

Interval-Sampling can be activated or deactivated through the Service-Menu.

### 4.2. Sterilization and Decontamination Before Use

If necessary, the rotor, protection cap, and air direction ring can be autoclaved at 121°C for 20 minutes. Alternatively, disinfection can be carried out using commercial disinfectants with a pH range between 5 and 8.

The housing of the RCS High Flow is made of polycarbonate (Makrolon™) and can be disinfected by spraying or wiping with commercial disinfectants which are harmless to polycarbonate. The suppliers of disinfectants can usually provide such information.


To protect the motor and the electronic circuits, exposure to excessive amounts of moisture, from immersion or over-use of sprays should be avoided. The instrument itself can not be autoclaved.

For the use of the RCS High Flow in aseptic environments it is recommended to use sterile sleeves (Art. No.940250). That will cover the instrument parts which are not autoclavable. Please take care that the outlet apertures of the air direction ring are not covered by the sterile sleeve.

### 4.3. Insertion of the Agar Strip

To insert the agar strip, first remove the protection cap by rotating it in a coun-

ter clockwise direction and then remove the rotor from the magnetic flange. Pull back the protective wrapper of the agar strip by approximately 4 cm

and remove the strip. Slide the strip into the rotor with the culture medium facing inwards. The top of the agar strip marked with the batch number and type should be flush with the guide plate of the rotor. Then replace the rotor onto the magnetic flange and replace the protection cap.

#### 4.4. Switching the instrument on

**Note:** For safety reasons the instrument can only be operated with the protection cap in place.

The instrument is switched on by pressing the [MENU/ON] key. By keeping the [MENU/ON] key pressed down at least 5 seconds, the version number of the instrument's internal software is displayed. The next display shows the serial number of the last rotor used.

If only one rotor is listed in the RCS High Flow rotor table, the volume selection is automatically activated. If more than one rotor is listed, a rotor can be chosen (see page 28).

#### 4.5. Setting the Sample Volume [☺]

The instrument can store a total of 10 sample volumes. The first seven memory locations are permanently set to the preset values of 10, 20, 50, 100, 200, 500 and 1000 liters in ascending order. In addition, three memory locations are available for entering individual volume settings between 1 and 1999 liters. It is recommended to choose a sample volume of 10 to 1000 liter.

The sample volume is adjusted using the [+] and [-] keys. A value one greater than

the value shown in the display is selected using the [+] key and one less than the value shown in the display using the [-] key.

Individually selectable volumes can be entered in three additional programmable

positions beyond the 1000 liter sample volume. For changing the individual volume refer to chapter 6 "Service-Menu".

The selection of a suitable sample volume depends on the expected concentration of microbes. Some examples are given in the following table:

Total count of microorganism expected in m <sup>3</sup> (CFU)	Pre-selected sample volume for RCS High Flow (Liters)	Number of expected colonies found on agar strip
0 to 10	1.000	0 to 10
10 to 200	500	5 to 100
200 to 500	200	40 to 100
500 to 1.000	100	50 to 100
1.000 to 2.500	50	50 to 125
2.500 to 10.000	20	50 to 200
10.000 to 20.000	10	100 to 200

#### 4.6. Starting the Instrument

After the required volume has been set, the air sampler is started by pressing

the [START/STOP] key on the keyboard of the RCS High Flow.

#### 4.7. Switching the Instrument Off

After the set volume of air has been sampled, the measuring process ends automatically. The value »0L« will appear in the display and an alarm will sound which can be switched off by pressing the [MENU/ON] key. The current measurement process can be interrupted if necessary by pressing the [START/STOP] key. To turn the instrument off, press and hold the [MENU/ON] key for 5 seconds. Alternatively, if left undisturbed the instrument will turn off automatically and the display will extinguish itself after eight minutes. It is advisable to check that sampling was carried out suc-

cessfully by switching the instrument on again using the [MENU/ON] key. If the sampling procedure was successful, the last used sampling volume will be shown in the display.

#### 4.8. Removal of the Agar Strip

To remove the agar strip, first remove the protection cap, then the rotor. With the help of the tab, pull the agar strip from the rotor and replace it into the protective wrapper with the culture media facing downwards. Seal the protective wrapper

using one of the supplied slide seals or adhesive tape to prevent the culture media from drying out during incubation.

Label the protective wrapper with the most relevant sample data to avoid mistakes of identification or use barcoding.

#### 4.9. Setting the Delay Time

If a Delay Time is required before starting the sampling procedure, this can be selected through the Service-Menu (see chapter 6.o).

#### 4.10. Power Supply

A fully charged 9.6 Volt battery allows the instrument to operate for approximately 1 hour. The recharging of the battery is also one hour, thus providing uninterrupted operation when using a spare battery. The battery is recharged with the battery recharger (Art. No. 940 370). The battery has to be placed into the compartment of the recharger labeled with 9.6 V and facing its »+« label towards the »+« label of the battery recharger.

**Note:** To avoid a reduction of battery capacity due to memory effects, it is recommended that the battery be recharged only when it is completely discharged or when the error message for low battery capacity is seen in the display.


The RCS High Flow Air Sampler can also be operated with a mains power supply (Art. No. 940 270: prim. 230/240 V, sec. 12 V / Art. No. 940 271: prim. 100/110 V, sec. 12 V) which can be plugged into the connection for mains adapter (No. 7, Section 1A.). For the use of this power supply a fully charged battery should also be placed in the RCS High Flow. The battery is required to start the sampling procedure. Using this combination the RCS High Flow can be used for several 100 measurements without the need for replacing the battery.

#### 4.11. Tripod and Table-Top Tripod

The RCS High Flow Air Sampler can be attached to a tripod which allows sampling at a height of approximately 3 m. For horizontal operation of the air sampler on very smooth surfaces it can be attached to a table-top tripod in order to minimize vibrations. Before the instrument is attached to the tripod adapter the protective screw located at the bottom of the RCS High Flow has to be removed.

#### 4.12. RS232 Interface / RCS Management Software

The RCS High Flow is equipped with an RS232 interface covered by a protective cap. It is recommended to remove the protective cap only when using the RS232. With help of the serial cable, which is delivered with the instrument, a transfer

Interface connection RS232

of instrument specific data to a personal computer is enabled. With the help of the RCS Management Software delivered with the instrument Interval-Sampling and quality safety levels can be selected and activated.

If activated, the QS module inhibits manual overwriting of the calibration factor on an instrument (QS level 2), the volume selection on an instrument (QS level 3) as well as the selection of one or more rotors by the user (QS level 4).

## 5. Evaluation of Results

The colonies or CFU are counted by direct visual inspection after incubation in the sealed plastic protective wrapper. The microbial count is to be stated with reference to the sampling volume.

The microbial count (CFUs) corresponding to the selected sample volume should be recorded as CFUs / sample volume. For CFUs per cubic meter (m<sup>3</sup>) the following calculation can be used:

$$\text{CFU/m}^3 = \frac{\text{CFU counted on agar strip}}{\text{sample volume (liter)}} \times 1000 \text{ (liter)}$$

For calculation of CFUs per cubic foot the following conversion factors apply: 1000 liter = 1 m<sup>3</sup> = 35.3 ft<sup>3</sup>; 1 ft<sup>3</sup> = 28.3 liter. To convert CFU/m<sup>3</sup> to CFU/ft<sup>3</sup> divide results by 35.

## 6. Service Menu

After switching the RCS High Flow Air Sampler on you are in the volume selection menu. Press the [MENU/ON] key to enter

the main Service Menu. To go back to the volume selection menu select the back arrow and press [MENU/ON] to confirm.


Service Menu

Press the [START/STOP] button to scroll through the Service Menu. In this menu you are able to select between: Bluetooth Data transfer, HYCON-ID, Delay Time, the setting of individual sampling volumes, Interval Sampling, and time setting.

### 6.1. Bluetooth Data transfer

After selecting the [Bluetooth] symbol for Bluetooth Data transfer you are able to enter the Bluetooth Data Menu. The Bluetooth Data transfer is used in conjunction with

Selection of Bluetooth Data transfer

Bluetooth Data transfer

HYCON-ID (Art. No. 940 360) and allows for an easy and secure means for transferring data to a personal computer. When the Bluetooth Data transfer is not in use the field under the disc symbol [OFF] must state „OFF“.


If it is displayed "ON", it must be changed into "OFF" by pressing the [+] or [-] keys.

Pressing the [MENU/ON] key returns back to the Service-Menu.

For activation of the wireless data transfer refer or other information regarding HYCON-ID refer to the operation manual.

### 6.2. HYCON-ID

After selecting the symbol representing HYCON-ID [PC] press the [MENU/ON] key to select. It is only necessary to select the HYCON-ID Service Menu if the Bluetooth Data transfer will be used.


Selection of HYCON-ID

HYCON-ID menu

Press the [MENU/ON] key to return to the Service Menu.

### 6.3. Delay Time

Choose the [Clock] to select Delay Time. Press the [START/STOP] key to enter the Delay Time in hours (either 0 or 1) and minutes. By


Selection of delay time

Delay Time menu

continuously pressing the [+] or [-] keys you are able to select and enter the desired Delay Time.

In order to activate the Delay Time change the field under the symbol by pressing the [+] or [-] key to switch between „ON“ and „OFF“.


If „ON“ is selected Delay Time is activated. Next, press the [MENU/ON] key to return back to the previous menu. Select the desired sampling volume and press [START/STOP].

Setting of Delay Time

Activation of Delay Time

Displayed in the upper right corner the Delay Time will begin counting down.

Volume display and Delay Time

### 6.4. Individual sampling volume

Select the [Volume] symbol and press [MENU/ON] to choose the required sampling volume. If necessary use the [+] and [-] keys to select one of the three open slots if choosing a sample volume other than one of the 7 pre-set volumes.

Selection of individual sample volume

Indiv. sample volume menu

*Selection of one of the individ. sampling volumes for setting*

Press the [START/STOP] and subsequently the [+] and [-] key to select the desired sampling volume.

*Setting of the individual sampling volume*

Press the [MENU/ON] key to return to the Service Menu.


In the volume selection menu the individual volume can be selected.

*Selection of an individual sampling volume in volume selection*

## 6.5. Interval-Sampling

The Interval Sampling feature on the RCS High Flow is used to extend a sample volume over a period of time. For example this feature is ideal for monitoring

*Selection of Interval-Sampling*

*Interval-Sampling menu*

a production run. This operation can only be set-up and activated through the RCS Management Software.

After establishing your parameters using the RCS Management Software, enter the Interval Sampling menu by selecting the [⏸] symbol and pressing the [MENU/ON] key on your RCS High Flow.


If activated it is possible to change between [ON] (active) and [OFF] (inactive) using the [+] and [-] keys.

Press the [START/STOP] key to activate. The display shows the number of intervals remaining in the upper left corner and the time between the intervals in the upper right corner.

After the unit has finished sampling, press the [MENU/ON] key to return to the Service Menu.

## 6.6. Setting Time

Select the [🕒] symbol to set the time.

*Selection of Time setting menu*

*Time setting menu*

By pressing the [START/STOP] - key you can navigate in the display between hours and minutes. Use the [+] and [-] - keys to highlight the selected field.

The date can only be changed using the RCS-Management Software. It is recommended to use the RCS Management Software setting both data and time.



Press the [MENU/ON] key to return back to the Service Menu.





*Selection time change*

*Changing time settings*

## 7. Display Messages/ Acoustic signals

Pos.:	Display	Symbol / Error number	Description
1.	Symbol flashes		Protection cup is missing
2.	Symbol flashes		Rotor is jammed
3.	Symbol flashes		Battery discharged
4.	Fixed		Malfunction of EEPROM
5.	Fixed		Malfunction: • Volume table / Delay-Time / SRAM – Data
6.	Fixed		Malfunction: • Calibration factor
7.	Fixed		Malfunction: • Rotor status
8.	Fixed		Malfunction: • Rotor expired
9.	Fixed		Malfunction: • QA-Level (Quality Control)
10.	Fixed		The RTC-second signal was not detected after start. Die RTC area check failed
11.	Fixed		During manual calibration, e.g. calculation check, a time conflict was detected with Interval-Sampling data
12.	Fixed		Malfunction during calculation check of Interval-Sampling data

Pos.:	Display	Symbol / Error number	Description
13.	Fixed		Write error
14.	Fixed		Write error
15.	Fixed		Expire date of checked rotor is reached
16.	Fixed		Malfunction of Bluetooth communication
17.	Fixed		Malfunction of Bluetooth Data transfer to PC
18.	Fixed		Malfunction during reading from data base of HYCON ID. Number of reading errors are displayed.
19.	Fixed		Malfunction during reading from data base of HYCON-ID. Number of errors are displayed.
20.	Fixed		HYCON ID data base is full.
21.	Message is shown for three seconds after start.		After start the actual data storage level of the HYCON-ID data base (>= 80%) is shown for three seconds.
22.	Fixed		This function is blocked (due to the selected by QA-level).

Various acoustic signals can sound which indicate either that the sampling process has been completed or that certain alarm conditions have been encountered:

 Sampling complete

 Alarm signal

## 8. Service and Maintenance

### 8.1. General - RESET


If part of the memory should be lost through problems arising from the power supply (e.g. removal of the battery during operation), this can be rectified by a General Reset. Memory loss becomes evident through the loss of one or more of the characters or individual segments on the display.

A General Reset is performed as follows:

With the instrument switched off press [MENU/ON] key and hold. After 10 seconds the menu of the General Reset is displayed.

Menu General Reset

Activation  
General Reset

Finished  
General Reset

By pressing the [+] or [-] keys the General Reset can be activated (display shows "ON"), press the [START/STOP] key to start the General Reset.

After finishing the General Reset the instrument displays "O.K."

Press the [MENU/ON] key to return back to the Service-Menu.

All permanently programmed and variable values for sampling volumes are reset to the factory values. The calibration correction factor entered into the instrument during calibration is not affected.

The selected time is not inactivated by a General Reset.

### 8.2. Calibration

The RCS High Flow stores a calibration expiration date for each rotor. The expiration date is automatically calculated when first switching on the instrument or during the performance of an electronic calibration. The expiration period of a calibration is 365 days. However, this period can be altered electronically using CaLibSo.

If a rotor was selected by pressing the [MENU/ON] key and the expiration date will soon expire, the serial number display of the rotor will blink. Confirm by press the [MENU/ON] key by which you can proceed to collecting the sample volume.

If a higher quality level (2-4) was selected with the RCS Management Software in QS management, the display will read an error message if a rotor with an exceeded calibration date was selected. By pressing the [MENU/ON] key you return to the rotor selection menu. By pressing the [+] and [-] keys another rotor can be selected. If no rotor with a valid expiration date is available, a calibration must be performed.

If the instrument is calibrated manually the last value e.g. 365 days is used as time period.

The Biotest RCS High Flow Air Sampler can be calibrated using the Biotest Anemometer (Art. No. 940 320) every 12 months in minimum from the date of first use.

Calibration is further recommended after dismantling the rotor.

**Note:** The procedure of calibration is described in the operating manual of the Anemometer. Due to the increased flow rate of the RCS High Flow an increased sample volume of 200 liters is recommended for measurement of the running time. In addition the formula described below for calculation of the correction factor should be followed.

If a discrepancy is found in the air volume sampled by the instrument, this can be corrected within in the range of ± 20% in increments of 2%.

The correction factor (Cf) to be set is derived from the following formula:

$$Cf = \frac{\text{selected sample volume (l)} \times 6000 \text{ (s/min} \times \%)}{\text{rotor flow rate (l/min)} \times \text{running time (s)}} - 100 \%$$

**Example:**

Result of differential measurement (= rotor flow rate)	<b>97 l/min</b>
Result of running time measurement for 200 l sample volume and a correction factor = 0 %	<b>122 s</b>
Selected sample volume	<b>200 l</b>


$$Cf = \frac{200 \times 6000}{97 \times 122} - 100\%$$

$$Cf = + 1,4\%$$

This means that the sampling time of the instrument must be increased by 2% or the correction factor set to ÷ 002 (= +2 %).

The setting of the correction factor is performed as follows:


After selection of the according symbol in the Service-Menu you enter the menu of the calibration factor by pressing the [MENU/ON] key.

Selection of Calibration factor menu

Calibration factor menu

By pressing the [START/STOP] key the calibration factor can be selected.

Selection of Calibration factor

Setting of Calibration factor

By highlighting the symbol [?] the calibration factor can be changed in double-steps between -20 and 20 using the [+] and [-] keys.

Press the [MENU/ON] key to return back to the Service-Menu.

**Note:** If several rotors are used with an instrument, the instrument must always be calibrated with each individual rotor. The correction factor determined for each rotor is stored in the memory of the air sampler and set up after the according rotor is selected.

### 8.2.1. Electronic calibration

Adding or deleting rotor numbers can be performed only electronically.

### 8.3. Cleaning of the Rotor

Occasional thorough cleaning of the rotor is recommended. Therefore the complete rotor can be cleaned in an ultrasonic water bath and rinsed with water afterwards or cleaned using commercially available disinfectant.

The rotor can also be dismantled into its individual parts after unscrewing the knurled nut and then cleaned carefully with a brush or in an ultrasonic waterbath. As this is a precision component which can influence the flow rate characteristics of the instrument, these parts must be handled with extreme care.

If several rotors are dismantled at the same time, it is of vital importance that the rotor housings and air direction cylinders are not interchanged, i.e., as rotor housings and air direction cylinders are numbered, they always bear identical numbers.

## 9. Accessories

The wide range of accessories available for the Biotest RCS High Flow Air Sampler makes the operation of the instrument even more convenient.

### 1. Agar Strips

#### 1.1. Agar Strips TC

(Art.No.941105); Casein Soybean Digest Agar for determination of total counts

#### 1.2. Agar Strips TSM

(Art.No.941135); modified Casein Soybean Digest Agar with neutralizers against disinfectants and growth supplements for determination of total counts of fastidious and sublethal damaged airborne microorganisms

#### 1.3. Agar Strips TC- $\gamma$

(Art.No.941115); double-wrapped,  $\gamma$ -irradiated Casein Soybean Digest Agar for determination of total counts in aseptic environments

#### 1.4. Agar Strips TCI- $\gamma$

(Art.No.941125); double-wrapped,  $\gamma$ -irradiated Casein Soybean Digest Agar with neutralizers for the determination of total counts in aseptic environments and also in peroxid-containing air

#### 1.5. Agar Strips YM

(Art.No.941196); Rose Bengal Agar for determination of yeasts and molds

#### 1.6. Agar Strips SDX

(Art.No.941205); Sabourand Dextrose Agar with modified Pharmacopoe formulation for determination of yeasts and molds

#### 1.7. Agar Strips SDX- $\gamma$

(Art.No.941215);  $\gamma$ -irradiated Sabourand Dextrose Agar with modified Pharmacopoe formulation for determination of yeasts and molds

#### 1.8. Agar Strips DG-18

(Art.No.941225); Dichloran Glycerole Agar (according to TRBA 430) for determination of yeasts and molds

#### 1.9. Agar Strips S

(Art.No.941405); Mannitol Salt Agar for determination of staphylococci

#### 1.10. Agar Strips C

(Art.No.941505); MacConkey Agar for determination of coliform bacteria

#### 1.11. Agar Strips PEN (Penase)

(Art.No.941705); Casein Soybean Digest Agar with Penase for determination of total counts in penicillin-production environments

#### 1.12. Agar Strips LAC (Lactamose)

(Art.No.941605); Casein Soybean Digest Agar with Lactamose for determination of total counts in antibiotic-production environments

### 1.13. Blank Strips

(Art.No.941605); sterile empty strips for self production of culture medium for special applications

The agar strips are produced under standardized conditions, thus enabling reproducible results to be obtained. Each package contains a batch specific quality control certificate.

### 2. Battery pack P 9.6

(Art.No.940275); 9.6 V spare rechargeable battery for uninterrupted operation of the RCS High Flow Air Sampler

### 3. Battery recharger

(Art.No.940370 - 230 V / 940371 - 110 V); for rapid charging of the battery pack P 9.6

### 4. Power supply

(Art.No.940270: prim. 230/240 V, sec. 12 V / Art.No.940271: prim. 100/110 V, sec. 12 V)

### 5. RCS Plus Anemometer

(Art.No.940320); portable air flow-rate measurement device for calibrating the Biotest RCS High Flow

### 6. Adhesive tapes

(Art.No.940325); for calibrating the RCS High Flow; adhesive tapes for complete sealing of the air inlet apertures on the rotor

### 7. Tripod

(Art.No.940330); for convenient operation of the RCS High Flow in heights up to 3 m

### 8. Table-top tripod

(Art.No.940335); for horizontal operation of the RCS High Flow on very smooth surfaces to minimize vibrations

### 9. Sterile sleeves

(Art.No.940250); for covering the RCS High Flow to protect the instrument against transmission contaminants

### 10. Air Outlet Ring

(Art.No.940244); spare part, autoclavable

### 11. Rotor, complete

(Art.No.940230); spare part, autoclavable

### 12. Protection cap

(Art.No.940415); spare part, autoclavable

### 13. Compressed gas adapter

(Art.No.940720); autoclavable adapter for monitoring microorganisms in compressed gases

### 14. HYCON - ID

(Art.No.940360); Set with barcodescanner for wireless data communication between RCS High Flow and PL.

## 10. Technical Data

**Sampling principle:** The Biotest RCS High Flow Air Sampler operates on the principle of impaction, whereby the air stream enters the rotor from the front of the instrument and the airborne microbes are separated onto the agar strip by centrifugal force. The air outlet passage is directed to the rear, parallel with the instrument

**Measurement range:** 1 -1999 liters (recommended measurement range 10 -1000 liters)

**Air flow-rate:** approx. 100 l/min

**Instrument precision:** ± 5%

**Sample volume:** 7 volumes in memory in ascending order of 10, 20, 50, 100, 200, 500 and 1000 Liters - 3 positions in memory, individually selectable from 1 to 1999 liters. The accuracy of sampling volume is achieved through calibration.

**Rotor speed:** approx. 8200 min<sup>-1</sup>

**Maximum allowable axial force on rotor shaft:** 30 Newtons

**Environmental conditions:** Temperature: 5 to 40°C; rel. humidity: 40 to 80%  
Power supply: rapid charging 9.6 V Ni/Cd battery (Art. No.940275); mains power supply (Art.No. 940270: prim. 230/240 V, sec. 12 V / Art.No.940271: prim. 100/110 V, sec. 12 V)

**Overvoltage Category:** II

**Pollution Degree:** 2

**Automatic power-off:** The instrument can be switched off by pressing the [ MENU/ON ] key in the volume setting display for five seconds or switches itself off automatically after 8 minutes.

#### Materials:

- **Housing:** polycarbonate; resistant to commercial disinfectants which are suitable for polycarbonate (Makrolon™)
- **Air direction ring:** polycarbonate (APEC®); autoclavable
- **Rotor:** anodized aluminium; autoclavable
- **Protection cap:** stainless steel; autoclavable

**Weight:** approx. 1500 g with battery pack (3.3 lbs)

#### Miscellaneous:

- Operation via keyboard panel with integrated display
- Selectable Delay Time
- Thread connection for tripod
- Can be calibrated using the Biotest Anemometer (traceable to PTB, Physikalisch Technische Bundesanstalt)
- Errors indicated via display and acoustic signal
- Memory function for residual sample volume and last sample volume processed
- Interface RS232, for connection with a personal computer via serial cable

As part of its continuous product improvement, Biotest reserves the right to amend technical specifications without prior notice.

## 11. Warranty

Biotest declares its instruments to be free from defects of materials and workmanship. The period of warranty is one (1) year commencing on the date of commissioning of the instrument as evidenced by means of an invoice or other document.

The warranty covers the rectification of all faults which occur as a proven direct result of faulty materials or workmanship during the warranty period.

The following items are not covered by the warranty:

1. Parts subject to normal wear and tear.
2. Malfunction caused by accidental or intentional damage, or the use of parts from other manufacturers, or through repairs or modifications carried out by persons not authorized by Biotest to carry out such work.
3. Defects which occur as a result of normal wear and tear or from using the equipment under conditions not specified in the instructions for use.

The user is not entitled to any remedy in law towards Biotest, howsoever arising, apart from those remedies mentioned here. Biotest accepts no other liability or warranty towards the user. The legal rights of the purchaser towards the seller remain unaffected by this.

**Biotest AG** • Landsteinerstr. 5  
63303 Dreieich • Germany

## 12. Waste Disposal

Expired instruments are electronic devices which should not be disposed with domestic waste. Biotest AG will deposit old instruments, when they are sent back carriage free.

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
## 1. Individual Parts and Functions

- ⑰ Battery compartment
- ⑱ Cover for battery compartment
- ⑲ Protective screw for stand mounting thread
- ⑳ Magnetic coupling
- ㉑ Bayonet lock for the protection cap
- ㉒ Knurled nut
- ㉓ Rotor housing
- ㉔ Rotor blades
- ㉕ Air direction cylinder
- ㉖ Magnet
- ㉗ Fan spindle
- ㉘ Interface RS 232

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Biotest Italia S.r.l.
Biotest (Schweiz) AG
Biotest Diagnostics Corp.