



FLOW CHEMISTRY



Beispielreaktionen:

+++ Grignaud +++ Alkylierungen +++ Wittig +++
Brominierungen +++ Chlorirung +++ Nitrierung
+++ Simmens-Smith +++ PMP +++ Aminierung +++
Mithylierung +++ Bouvault Aldehyd-Synthese
+++ Mitsunobu-Reaktion +++ Suzuki-Coupling +++
+++ NaBh₄-Reduktion +++ Sonogashira-Coupling +++

Allgemeines


Mikroreaktionstechnik zählt gegenwärtig zu einer der innovativsten Technologien auf dem Gebiet der chemischen Synthese und verwandter Gebiete. Mikrostrukturierte Reaktoren halten jetzt Einzug in die Labore. Weltweite Forschungsaktivitäten zeigen eindrucksvoll, dass sich durch den Einsatz von Mikroreaktoren eine Vielzahl von reaktionstechnischen Vorteilen in chemischen Prozessen und Reaktionen erzielen lassen. YMC Europe bietet eine vollständige Produktreihe, größtenteils automatisierter Mikroreaktoren für Forschung und Produktion an. Mit den Produktserien KeyChem und CYTOS kann der gesamte Syntheseworkflow vom Forschungslabor bis zum Technikumsmaßstab wirkungsvoll unterstützt werden.

Ein umfangreiches Baukastensystem für die ersten Schritte im Forschungslabor bis hin zur Produktion unterstützt sowohl Einsteiger, als auch erfahrene Anwender bei der Umsetzung neuer Synthesestrategien. So sind erste Orientierungsexperimente mit dem KeyChem Basic System einfach und schnell umsetzbar. Mit den KeyChem-L Modulen stehen automatisierte, softwaregesteuerte Systeme mit umfangreichem Datenlogging sowie einer breiten Auswahl an Reaktoren und Zubehör zur Verfügung. Die CYTOS 200 und CYTOS 2000 Serie stellt den Einstieg in den Technikums- und Produktionsmaßstab dar.

Die Auswahl der Reaktoren erfolgt aufgrund der Erfordernisse der durchzuführenden Reaktionen und aufgrund der Art und Weise, wie der Kunde die Anlage einsetzen möchte. Die Fördermodule sind ebenfalls von den eingesetzten Reaktoren abhängig. Im Bereich der flüssigen Medien erstreckt sich das Spektrum der eingesetzten Pumpen von analytischen HPLC Pum-

pen bis hin zu hochdruckresistenten Präzisions-spritzenpumpen. So ist es mit Präzisions-spritzenpumpen möglich, konstante Stoffströme bis in den Nanoliterbereich kontinuierlich auch gegen hohe Drücke zu fördern. Die Module dieses Systems arbeiten auch unter Druck stabil und störungsfrei, selbst bei der Förderung viskoser Medien. Unterschiedliche Mischer und Weir modules von wenigen Mikrolitern bis hin zum großvolumigen Milliliter-Modul sind mit unterschiedlicher Geometrie verfügbar. Innovative Peltiertechnik ermöglicht bei den KeyChem-Systemen Heizen und Kühlen ohne Kryostat oder zusätzliche Wärmeträgermedien.

YMC Mikroreaktoren leisten einen Beitrag zum Konzept der „Green Chemistry“, indem sie Anwendern ermöglichen, ressourcenschonende Synthese- und Produktionsverfahren zu finden und in den Produktionsprozess zu integrieren.



SYSTEMS

KeyChem-L:

Das System KeyChem-L bietet Mikroreaktionstechnik mit kompaktem Design, einer Spritzenpumpe und Temperaturregelung durch Peltier Elemente.

Produktmerkmale

- Ideale Positionierung im Labor durch modulare Bauweise; PC, Pumpe und Thermoelement sind voneinander getrennt.
- Das System kann auf die individuellen Ansprüche jedes Kunden angepasst werden
- Benutzerfreundliche Software zur Steuerung
- Hervorragende Durchmischung und Temperierung
- Und nicht zuletzt, ein günstiger Preis!



KeyChem Lumino:

Für photochemische Reaktionen können KeyChem-Systeme einfach und schnell aufgerüstet werden:

KeyChem-Lumino ist ein Photomikroreaktor, wahlweise aus rostfreiem Stahl (SUS316) mit Deckel aus Quarzglas oder TiO₂-beschichtet, der von einer abgeschirmten Hochleistungslichtquelle bestrahlt wird.

KeyChem-Lumino Photomikroreaktoren sind mit UV-Quellen verschiedener Wellenlängen und Energieeinträgen lieferbar:

Niederdruck-Quecksilber-Lampe mit 254 nm Wellenlänge

Excimer-Lampe 308 nm oder


regelbare 6-fach LED's mit 365 nm.



KeyChem-L/HP:

KeyChem-L-HP ist ein kontinuierlich arbeitendes Autoklavensystem für alle Reaktionen, die unter Druck stattfinden können. Im Mitteldruckbereich sind 30 bis 50 bar möglich, Hochdruck-Reaktoren auf Anfrage lieferbar. Alle Reaktoren und Verweilschleifen sind in Edelstahl oder Hastelloy verfügbar.

Die Temperierung erfolgt über peltiergesteuerte Thermostaten, die schnelles Heizen und Kühlen ermöglichen.



CYTOS:

CYTOS ist ein modulares System zur Produktion größerer Mengen, mit einer Größenordnung von einigen 100 Gramm bis Kilogramm – je nach Syntheseverfahren und Auslegung der Anlage. Durch Verschalten mehrerer CYTOS Module (Numbering – up) sind größere Stoffumsätze erreichbar. Durch Parallelisieren kann aus einem einzelnen CYTOS Reaktor ein leistungsfähiges Produktionssystem aufgebaut werden. Dieses Konzept hat den Vorteil der größten möglichen Flexibilität im Anlagenbetrieb. Fällt ein Reaktor aus, kann ein Reservesystem dessen Funktion übernehmen, ohne dass die Funktionalität der Anlage beeinträchtigt wird. Für CYTOS steht eine abgestimmte Auswahl an Dosier- und Temperier-Komponenten zur Verfügung, die über eine intuitiv zu bedienende Software mit Datenaufzeichnung gesteuert werden können.



FLOW CHEMISTRY

FRX – Flow Chemistry



Die FRX-Reihe von Syrris macht FlowChemistry kosteneffektiv, einfach in der Handhabung, verständlich und gleichzeitig kostengünstig. Die FRX-Reihe besteht aus mehreren einzelnen Modellen in modularer Bauweise, so dass diese unkompliziert zu einfachen aber leistungsstarken FlowChemistry-Systemen zusammengestellt werden können.

Africa – automatisierte Flow Chemistry



Africa ist das komplexeste und flexibelste Mikroreaktorsystem, welches verfügbar ist. Reaktionsoptimierungen, programmierte Synthesen oder präparative Synthesen können auf dem Africa-System mit anschließender Aufbereitung und Analyse vollautomatisch realisiert werden.

YMC Auftragssynthese und Entwicklung

Das Labor in Dinslaken (Nähe Düsseldorf) unterstützt Anwender wirkungsvoll in der Methodenentwicklung und -optimierung. Dort werden auch Auftrags-synthesen, (präparative) Auftrags-trennungen mittels HPLC, Schulungen und Seminare durchgeführt.



YMC Europe GmbH

Schöttmannshof 19
D-46539 Dinslaken
Germany
TEL. +49(0)2064 427-0, FAX +49(0)2064 427-222
www.ymc.de

YMC Co., Ltd.

YMC Karasuma-Gojo Bld. 284 Daigo-cho,
Karasuma Nisiiru Gojo-dori Shimogyo-ku,
Kyoto 600-8106 Japan
TEL. +81(0)75-342-4515, FAX +81(0)75-342-4550
www.ymc.co.jp

YMC America, Inc.

941 Marcon Boulevard Suite 301
Allentown, PA18109 USA
TEL. +1-610-266-8650, FAX +1-610-266-8652
www.ymc-america.com

YMC India Ltd.

CX - 07, 3rd Floor, Lobe - 1,
Tower - A, The Corenthum, Plot No- A-41,
Sector - 62, Noida - 201301 (UP) India.
TEL. +91(0)120-4276020 - 25, FAX +91(0)120-4276026
www.ymcindia.com

In Österreich:


OmniChrom HandelsgmbH

Frömmelgasse 40-42/15 · A-1210 Wien
TEL. +43(0)1-2785701, Fax +43(0)1-2727574
www.ymc.de

In der Schweiz:

STAGROMA

Chr.-Merian-Ring 31a · CH-4153 Reinach
TEL. +41(0)61-7178717, Fax +41(0)61-7178710
www.stagroma.com



Chemical **Innovation**

Keyboard Chemistry®

Version 10.5

Microreactor general catalog

YMC Co., Ltd.

Microreactor Systems
Reactor-Residence time units
Thermostats
Pumps
Accessories-Consumable parts

YMC aims to create “Keyboard Chemistry[®]”.

“Keyboard chemistry[®]” has been coined by YMC by combining the words “Keyboard” (as in computer control) and “Chemistry”, to represent the new concept of a chemical reaction methodology that YMC is creating.

Central to this concept is the microreactor system, which allows computer-controlled flow reaction and enables the setting of all the conditions of a chemical reaction from start to finish, via the keyboard of a computer. Many innovative reaction parameters, which can not be achieved when using conventional flasks and reaction vessels now become possible.

By extending the concept of Keyboard Chemistry[®] from chemical reactions to peripheral technologies such as pre-treatment, post-reaction analysis, separation and purification, YMC aims to change the landscape of chemistry forever.


YMC supplies the microreactor equipment which brings this innovation for R & D and

What is a microreactor ?

The term “microreactor” is the generic name given to devices that enable continuous flow reactions which micro-fabricated modules (mixer, heat exchanger, etc.) with typical thicknesses below 1 mm.




YMC is ready to develop apparatus which satisfies customers’ needs and provides the optimized reaction conditions from laboratory to plant scale.





Company Information

Company name	YMC Co., Ltd.
Founded	January 28th, 1980
Activities	Production and sales of packing materials and packed columns for high performance liquid chromatography (HPLC), Custom purification and synthesis, Manufacture and sales of microreactor systems and devices
Head office address	YMC Karasuma-Gojo Bldg., 284 Daigo-cho, Karasuma Nishiiru Gojo-dori, Shimogyo-ku, Kyoto 600-8106, Japan
CEO	Ryuji YAMAMURA
Capital	627 million yen
Annual turnover (consolidated)	4254 million yen (2010)
No. of employees (consolidated)	246 (as of Mar. 2010)
Locations	YMC Co., Ltd. Komatsu Works (Komatsu City, Ishikawa prefecture) YMC Co., Ltd. Tokyo Sales Office (Chiyoda-ku, Tokyo) YMC Microreactor Research Center (Kumiyama-cho, Kyoto.) Shanghai representative office
Affiliated companies	Aishin Kohki Co., Ltd. YMC Europe GmbH YMC America Inc. YMC India (Pvt.) Ltd.



Head Office.



Microreactor Research Center



Tokyo Sales Office.



Komatsu Works.

History of Microreactor Business


Apr. 2005	Become an exclusive distributor of CPC (Germany) microreactors for Japan
Aug. 2007	Establish YMC Microreactor Research Center
Oct. 2007	Purchased all assets of CPC (Germany)
Nov. 2007	Became an exclusive distributor of CORNING SAS Advanced Flow Reactor for Japan
Apr. 2008	Start sales of CPC CYTOS [®] series microreactors made in Japan Start sales of the YMC original microreactor (KeyChem [®] -L)
Aug. 2008	Investment in Syrris, a UK microreactor manufacturer



KeyChem[®]-L



CYTOS[®]



Corning product



Syrris product

CONTENTS

Microreactor Systems

A

- 1 KeyChem® -L
- 2 CYTOS® -200
- 3 CYTOS® -2000
- 4 KeyChem® -Lumino
- 5 KeyChem® -H
- 6 KeyChem® -Basic
- 7 Microreactor simple kit

Reactors • Residence time units

B

- 1 KeyChem® Mixers Hadar (Y type •SUS/Hastelloy/Glass/PTFE)
- 2 KeyChem® Mixers Deneb (Helix type •SUS/Hastelloy/Glass)
- 3 KeyChem® Mixers Spica (Static type •SUS/Hastelloy/Glass)
- 4 KeyChem® Mixer Rigel (Y type RTU 4.5 mL •SUS)
- 5 KeyChem® Residence time units Sirius (1.5 mL/3.0 mL/4.5 mL •SUS)
- 6 CYTOS® -200 Reactors Sargas (Y type •SUS/Glass)
- 7 CYTOS® -200 Reactors Regulus (CYTOS® -M type •SUS)
- 8 CYTOS® -200 Residence time units Laminar&Capillary RTU (SUS/Hastelloy)
- 9 CYTOS® -2000 Reactors (CYTOS® -L type •SUS/Hastelloy)
- 10 CYTOS® -2000 Residence time units Laminar&Capillary RTU (SUS/Hastelloy)
- 11 KeyChem® -Lumino Reactor Polaris (SUS)
- 12 KeyChem® -Lumino Reactors Polaris (TiO₂)

Thermostatic units

C

- 1 Coupling type Thermo Stages (Low temperature version • 2 type/3 type/4 type)
- 2 Coupling type Thermo Stages (High temperature version • 2 type/3 type/4 type)
- 3 Stand-alone type Thermo Stage (Low temperature version • Mixer/RTU)
- 4 Stand-alone type Thermo Stage (High temperature version • Mixer/RTU)
- 5 KeyChem[®]-Lumino Thermo Stage

Pump

D

- 1 MR2 Syringe pump (1-5 way valvex4)
- 2 MR2 Syringe pump (1-5 way valvex2)、(1-2 way valvex2)

Accessories • Consumable parts

E

- 1 Tubing connector 1/16 inch (pack of 10) Stainless steel
- 2 Tubing connector 1/16 inch PEEK • M6HN
- 3 Luer lock Connector
- 4 KeyChem[®]-L Tubing connector (Low temperature version)
- 5 KeyChem[®]-L Tubing connector (High temperature version)
- 6 Connector 1/16 inch (pack of 10) ETFE
- 7 EASYFIT10φ(pack of 10)
- 8 Housing connector (Stainless steel tubing version)
- 9 Syringe adapter (female type)
- 10 Luer adapter (male type)
- 11 Ferrule 1/16 inch (pack of 10) PTFE
- 12 Ferrule 1/16 inch (pack of 10) Stainless steel
- 13 Ferrule 1/16 inch (pack of 10) PTFE • GF
- 14 Back pressure regulator 5 psi
- 15 Back pressure regulator 20 psi
- 16 Back pressure regulator 40 psi

- 17 In-line check valve
- 18 In-line filter set (includes 10 filters)
- 19 Suction filter (PP)
- 20 Suction filter (SUS)
- 21 T type reactor·inside diameter 0.5 mm (Stainless steel)
- 22 T type reactor·inside diameter 1.0 mm (Stainless steel)
- 23 T type reactor·inside diameter 0.25 mm (PEEK)
- 24 T type reactor·inside diameter 0.5 mm (PEEK)
- 25 T type reactor·inside diameter 1.0 mm (PEEK)
- 26 PTFE tubing (Outside diameter 1/16 inch·inside diameter 0.25 mm·length 10 m)
- 27 PTFE tubing (Outside diameter 1/16 inch·inside diameter 0.5 mm·length 10 m)
- 28 PTFE tubing (Outside diameter 1/16 inch·inside diameter 0.75 mm·length 10 m)
- 29 PTFE tubing (Outside diameter 1/16 inch·inside diameter 1.0 mm·length 10 m)
- 30 Stainless steel tubing (Outside diameter 1/16 inch·inside diameter 0.5 mm·length 10 m)
- 31 Stainless steel tubing (Outside diameter 1/16 inch·inside diameter 1.0 mm·length 10 m)
- 32 KeyChem[®] Mixer housing
- 33 KeyChem[®] Mixer housing (Acrylic)
- 34 KeyChem[®] Residence time unit housing
- 35 CYTOS[®]-200 Reactor housing
- 36 CYTOS[®]-200 Clamping unit
- 37 KeyChem[®]-Lumino Reactor housing
- 38 Hamilton Syringe attachment kit for MR2 pump (2.5 mL/1 mL)
- 39 Hamilton Syringe attachment kit for MR2 pump (5 mL)
- 40 Hamilton Syringe attachment kit for MR2 pump (10 mL)
- 41 Hamilton Syringe attachment kit for MR2 pump (25 mL)
- 42 Disposable Syringe attachment kit for MR2 pump (2 mL)
- 43 Disposable Syringe attachment kit for MR2 pump (5 mL)
- 44 Disposable Syringe attachment kit for MR2 pump (10 mL)
- 45 Disposable Syringe attachment kit for MR2 pump (20 mL)
- 46 KLOEHN Syringe attachment kit for MR2 pump
- 47 Reservoir base for MR2 pump
- 48 Reservoir base with pressure meter for MR2 pump
- 49 O-ring for CYTOS[®] (Heat fluid line)
- 50 O-ring for CYTOS[®] (Reaction line)
- 51 STACKING DRY BATH
- 52 UV light source unit for KeyChem[®]-Lumino (365 nm)
- 53 Low pressure mercury lamp for KeyChem[®]-Lumino (254 nm)
- 54 Excimer lamp for KeyChem[®]-Lumino (308 nm)



- KeyChem[®]-L
- CYTOS[®]-200
- CYTOS[®]-2000
- KeyChem[®]-Lumino
- KeyChem[®]-H
- KeyChem[®]-Basic
- Microreactor simple kit

A

Microreactor Systems



Description

Flow chemistry apparatus of compact design with a syringe pump and temperature control by a Peltier thermostat.

Features

- The separate construction of each component (PC, syringe pump and thermostat) allows optimum positioning.
- It allows system configuration for the user's needs.
- Compact design allows for easy installation.
- User-friendly software
- Excellent mixing and temperature control
- Affordable pricing



Specification

	Throughput	Several hundred g/day
Pump	Type	Syringe pump
	Syringe size	1 mL - 25 mL
	Flow rate range	0.7 µL/min - 55 mL/min (recommended flow rate range : 10 µL/min - 10 mL/min)
Temperature unit	Temp. type	Peltier control
	Range	Approx. -15 - 80 °C (Standard type) • High temperature type is available.(ambient temperature – 150 °C)
Reactor unit	Mixer type	Y type, Static type, Helix type
	Residence time unit	Volume : 1.5 mL, 3.0 mL, 4.5 mL
	Reactor material	Stainless steel 316, Hastelloy, Glass, PTFE

Mixers

Three types of micromixer and three kinds of Residence time units are available as standard.



Thermostatic units

Thermo stages controlled by a Peltier unit covering a wide range of temperature are available to suit the user's needs.



Fluid pumps

Three different brands of syringes are used as standard with easy mounting/removal. Continuous flow can be achieved by using two syringes with alternate strokes (optional).

Description

Optimized flow chemistry apparatus for R&D use with a syringe pump and an independent thermostat with a wider range of temperature control

Features


- Precise temperature control by excellent heat exchange efficiency
- It is an improved version of CPC CYTOS[®]-M.
- Robust design of long term performance
- Full PC control

Specification

	Throughput	Several hundred g/day
Pump	Type	Syringe pump (Piston pump is also available.)
	Syringe size	1 mL - 25 mL Syringe
	Flow rate range	0.7 μ L/min - 55 mL/min (recommended flow rate range : 10 μ L/min - 10 mL/min)
Temperature unit	Temp. type	Heat fluid circulation (A thermostat is not included as standard supply.)
	Range	-20 - 170 $^{\circ}$ C (Depending on the thermostat used)
Reactor unit	Mixer type	Sargas(Y type), Regulus(CYTOS [®] -M type)
	Residence time unit	Volume : 0.9 mL, 1.8 mL, 4.5 mL
	Reactor material	Stainless steel 316, Glass

PC control

User friendly software is provided to set, and monitors the pressure and the temperature. It also allows saving and analyzing of the experimental data.

Fluid pumps

Three different brands of syringes are used as standard with easy mounting/removal. Continuous flow can be achieved by using two syringes with alternate strokes (optional).

Description

Versatile flow chemistry apparatus covering a range from laboratory use to production scale with high performance piston pumps(0.1 mL - 20 mL/min)

Features

- Optimized apparatus for process development from research use to bench scale synthesis
- It is an improved version of CPC CYTOS®-L.
- Easy optimization based on accurate parameters setting
- Multistep reaction is available due to external modules.
- Continuous flow is available from piston pumps.
- Reactor volume is 10 times larger than that of CYTOS®-200 allowing easy scale-up to be achieved.




Specification

	Throughput	- Several kg/day
Pump	Type	Piston pump
Temperature unit	Temp. type	Heat fluid circulation (A thermostat is not included as standard supply.)
	Range	-20 - 170 °C (Depending on the thermostat used)
Reactor unit	Reactor	CYTOS®-L type
	Residence time unit	Volume : 15 mL, 60 mL
	Reactor material	Stainless steel 316, Hastelloy

Dedicated software

User friendly software is provided by way of the touch panel and PC control with a simple display. It allows monitoring of the temperature and the pressure and allows saving and analyzing of the experimental data.



Large volume microreactor

The reactor has the same design configurations as the CYTOS®-200's with the equivalent heat exchange and mixing efficiency. This allows easy scale-up with the same reaction conditions to be achieved. Several sizes of residence time units are available for different reaction



4 Microreactor System KeyChem[®] series KeyChem[®] -Lumino

Model number MicroreactorSystem-LU

Description

Flow chemistry apparatus dedicated for photo-reactions

This system can connect every pump (Please see the pumps section).

Features

- Excellent light irradiation efficiency
- Various light sources
(LED light source, low pressure mercury lamp, excimer lamp)
- Precise temperature control possible
- The reactor can be disassembled
- Compact design for optimization of work space



Specification


Light source	
Lamp (6 LEDs in light emitting unit)	
Light type	Ultraviolet emission LEDs
Peak wavelength	365 nm (25 °C)
Life time	500 hours in continuous operation
Light power	Max. 250 mW (each)
Power requirements	AC 100-240 V, 100 W
Dimensions	(WxDxH)
	Light emitting unit : 55x57x40 mm Power source unit : 50x250x161 mm
Reactor unit	
material	Stainless steel 316, Quartz glass(cover)
Channel width	1.0 mm
Channel depth	0.2 mm
Channel length	560 mm
Volume	110 μL
Dimension	φ50×20 mm (excluding projection)
Temperature unit	
Power requirements	AC 100 - 240 V, 100 W
Temp. type	Peltier unit
Range	-15 - 80 °C
Precision	±1 °C
Dimension	(WxDxH) 120x230x205 mm
Weight	Approx. 4 kg

High performance LED light source


Uses a high performance LED as light source
The power output can be adjusted
Light dose not leak outside reactor.



LED light source (365 nm)




Excimer lamp (308 nm)



Light sources

Three light sources available: LED light source, low pressure mercury lamp, excimer lamp. For additional information, please see the accessories section.

Low pressure mercury lamp (254 nm)



Reactor

A reactor which is coated on surface with TiO₂ as a light catalyst is also available (not the stainless steel reactor).



Description

High efficiency and safe hydrogenation is now a reality. The hydrogen and the substrate are mixed in a gas-liquid micromixer, and flow into the reactor column filled with catalyst.

Features

- New state-of-the-art hydrogen supplying system
(Using metal hydride)
- Catalyst-filled column (Pd polymer bound)
- Various catalysts can be selected
- Column is available in a wide range of sizes.

Specification

Temp. Range (column oven)	ambient temperature - 100 °C
Pressure resistance	0.98 MPa
Back pressure (Max)	0.7 MPa
Flow rate range	0.1 mL/min - 10 mL/min
Power requirements	AC 100 - 240 V, 50/60 Hz, 300 W
Dimensions	(WxDxH) 380x335x375 mm including the display 380x335x530 mm



For added flexibility, options include expanded ranges of flow, temperature, and pressure.

Please enquire.

*System option include expanded ranges for temperature, pressure, and flow rate Please inquire.

Hydrogen absorbing alloy cylinder

Hydrogen storage on metal hydride is not at high-pressure, so it can be used safely even inside laboratories.

Reduces running costs



Hydrogen absorbing alloy cylinder (Standard size)

Dimension	(O.D.xH) 54x270 mm
Hydrogen storage capacity	Approx. 220 L
Weight	Approx. 2.1 kg
Materials	aluminum alloy

Cylinder is available in various other sizes

Reaction column catalyst

The reaction column is filled with a dedicated hydrogenation polymer-supported Palladium catalyst. The catalyst is reusable.

Catalyst-filled column

Filling catalyst	Pd polymer bound, Pd, Pd/C
Standard column size	(I.D.xL) φ4.6x30 - φ10.0x150 mm

Columns can be supplied with other catalysts; and can also be supplied in sizes other than the standard size. Please enquire.




Description

Entry level unit for flow reaction based on KeyChem®-L
 Only manual operation is available.

Features

- Simple manual control
- Compact design
- Oil-free temperature control for safety operation
- Affordable pricing



Specification

Pump	Type	Syringe pump
	Syringe size	1 mL - 25 mL syringe
	Flow rate range	0.7 µL/min - 17 mL/min
Temperature unit	Temp. type	Peltier unit
	Range	-15 - 80 °C (Standard) High temperature type is available. (ambient temperature - 150 °C)
Reactor unit	Mixer	Y type, Static type, Helix type
	Mixer material	Stainless steel 316, Hastelloy, Glass, PTFE

Benefits of KeyChem®-Basic

- Compact space
- Simple reaction parameter setting
- Low price
- Commercially available pumps are usable.



Thermostatic unit















Independent control
 Temperature range : -15 - 80 °C or
 ambient temperature - 150 °C

Fluid pumps

Two types of syringe can be used.
 Flow rate range : 0.7 µL/min - 17 mL/min
 Syringe size : 1 - 25 mL



In addition to the standard specifications, the following units are available as options depending on users' requirements.

Pumps		Mixers • Residence time units		Thermostatic units	
MR2 syringe pump		KeyChem [®] Mixer Hadar (Y type)		Independent Thermo Stage (Low temperature version, Mixer)	
	Pump with four syringes Suction/Discharge (independently movable) Usage size/1 - 25 mL PC control		Material:Stainless steel 316 Channel structure:Y type		Range:-15 - 80 °C Oil-free thermostat using Peltier control for mixer
		KeyChem [®] Mixer Deneb (Helix type), Spica (Static type)		Independent Thermo Stage (Low temperature version, RTU)	
			Material:Stainless steel 316 Channel Structure: Deneb:3 dimensional Spica:repeating splitting and joining		Range:-15 - 80 °C Oil-free thermostat using Peltier control for residence time unit
		KeyChem [®] Mixer Glass		Independent Thermo Stage (High temperature version, Mixer)	
			Material:Glass Channel structure: Y type, Helix type, Static type		Range: ambient temperature - 150 °C Oil-free Thermostat by heater control for mixer
		KeyChem [®] Mixer Hastelloy		Independent Thermo Stage (High temperature version, RTU)	
			Material:hastelloy Channel structure: Y type, Helix type, Static type		Range: ambient temperature - 150 °C Oil-free Thermostat by heater control for residence time unit
		KeyChem [®] Mixer Hadar (Y type) PTFE		Housing kit (KeyChem [®] Mixer)	
			Material:PTFE Channel structure:Y type *Possible to disassemble		Housing of KeyChem [®] mixer for dipping
		KeyChem [®] Mixer Rigel (Y + 4.5 mL type)		Housing kit (KeyChem [®] RTU)	
			Material:Stainless steel 316 Channel structure:Y type + 4.5 mL residence time unit		Housing of KeyChem [®] RTU for dipping
		KeyChem [®] RTU Sirius (1.5, 3.0, 4.5 mL)			
			Material:Stainless steel 316 Plate of residence time unit for taking reaction time		

Yellow area : KeyChem[®]-Basic Standard set (at last page)



Description

Simple microreactor kit for beginners
(exclude the pump)

Features


- Allows investigation of microreactions with minimum equipment




Specification

Mixer unit	Mixer	Y type, Static type, Helix type
	Mixer material	Stainless steel 316, Hastelloy, Glass, PTFE

Parts detail



	Product	Model number
1	KeyChem® Mixer Hadar (Y type)	KC-M-Y-SUS
2	KeyChem® Mixer Spica (Static type)	KC-M-S-SUS
3	KeyChem® Mixer Deneb (Helix type)	KC-M-H-SUS
4	Connector 1/16 inch PEEK·M6HN (3 connectors)	YMC-P-0002
5	Connector 1/16 inch (pack of 10) ETFE	YMC-P-0006
6	Syringe adapter (female type) (2 adapters)	YMC-P-0007
7	Ferrule 1/16 inch (pack of 10) PTFE·GT	YMC-P-0012
8	PTFE tubing (O.D.:1/16 inch, I.D.:0.5 mm, Length:10 m)	YMC-P-0025
9	KeyChem® mixer housing	YMC-P-0030-01



- KeyChem® Mixer•Residence time unit
- CYTOS-200® Reactor•Residence time unit
- CYTOS-2000® Reactor•Residence time unit
- KeyChem®-Lumino Reactor

B

Reactors - Residence time units



1 KeyChem® series Mixers Hadar (Y type)

Model numbers see below


Description

Simple Y-shaped channel design for entry-level microreaction research

Features

- Excellent heat exchange efficiency
- Suitable for entry-level to microreaction systems
- Adaptable for particle synthesis
- Suitable for reactions with acidic liquids (Hastelloy, Glass, PTFE)

Specification



Size LengthxWidthxThickness(mm)	Volume(μL)		Channel size Min. WidthxDepth(mm)	Material	Model number
	Preheating	Mixing			
30x30x1.3	17	1.4	0.5x0.1	SUS316	KC-M-Y-SUS
30x30x1.8	34	2.8	0.5x0.2	Hastelloy	KC-M-Y-HC
30x30x3.2	14	1.2	0.5x0.1	Glass	KC-M-Y-G
30x30x2.5	12	3.8	0.5x0.5	PTFE	KC-M-Y-P

PTFE : different channel design



2 KeyChem® series Mixers Deneb (Helix type)

Model numbers see below


Description

Stirring mixer which has three dimensional spiral-like channel design

Features

- Suitable for organic synthesis
- Excellent heat exchange efficiency
- Suitable for reactions of acidic liquids (Hastelloy, Glass)

Specification



Size LengthxWidthxThickness(mm)	Volume(μL)		Channel size Min. WidthxDepth(mm)	Material	Model number
	Preheating	Mixing			
30x30x1.8	57	32	0.2x0.2	SUS316	KC-M-H-SUS
30x30x1.8	39	22	0.2x0.2	Hastelloy	KC-M-H-HC
30x30x3.2	52	33	0.2x0.2	Glass	KC-M-H-G

PTFE : Please enquire.



3 KeyChem® series Mixers

Spica (Static type)


Model numbers see below

Description

Laminar flow mixer with a series of division and combination of flow occurs

Features

- Excellent mixing efficiency
- Suitable for reactions of acidic liquids (Hastelloy, Glass)



Specification

Size LengthxWidthxThickness(mm)	Volume(μL)		Channel size Min. WidthxDepth(mm)	Material	Model number
	Preheating	Mixing			
30x30x1.8	57	20	0.2x0.2	SUS316	KC-M-S-SUS
30x30x1.8	39	15	0.2x0.2	Hastelloy	KC-M-S-HC
30x30x3.2	56	16	0.2x0.2	Glass	KC-M-S-G

PTFE : Please enquire.



4 KeyChem® series Mixer

Rigel (Y type RTU·SUS)

Model number KC-MR4.5-Y-SUS

Description

A reactor with a Y type mixer and a residence time unit (RTU) of 4.5 mL for organic reactions requiring longer reaction times

Features

- Compact design with combined mixer and RTU
- Requires only one Thermo Stage



Specification

Size LengthxWidthxThickness(mm)	Volume (μL)			Channel size Min. WidthxDepth(mm)
	Preheating	Mixing	Residence	
50x66x6.5	42	2.7	4,600	0.5x0.1



5 KeyChem® series Residence time units

Sirius (SUS)

Model numbers see below

Description

Plate type RTU for longer reaction times.

A dedicated Thermo Stage is needed.

Features

- Three different volumes
- Excellent heat exchange efficiency due to stainless steel 316 construction



Photo : Sirius15

Specification

Product name	Volume (mL)	Size LengthxWidthxThickness(mm)	Channel size WidthxDepth(mm)	Model number
Sirius15	1.5	50x66x3.7	2x0.6	R006-50-S6
Sirius30	3.0	50x66x5.1	2x1.3	R005-50-S6
Sirius45	4.5	50x66x6.5	2x2	R004-50-S6



6 CYTOS®-200 series Reactors

Sargas (Y type)

Model numbers see below

Description

Y type reactor for CYTOS®-200 with internal heating zones.

Stainless steel unit has internal heat fluid circulation zones for wider temperature control range.

Glass unit is suitable for handling acidic liquids.

Features

- Simple Y-shape mixing
- Excellent heat exchange efficiency
- Wider range of temperature control by heat fluid circulation



SUS316



Glass

Specification

Size LengthxWidthxThickness(mm)	Volume (μL)		Channel size Min. WidthxDepth(mm)	Material	Model number
	Preheating	Mixing			
50x130x6.4	94	42	0.5x0.2	SUS316	C200-Y-SUS
50x130x6.2	170	73	0.5x0.4	Glass	C200-Y-G

7 CYTOS[®]-200 series Reactor

Regulus (CYTOS[®]-M type · SUS) Model number C200-M-SUS

Description

Original CPC design for CYTOS[®]-200 reactor with internal heat fluid circulation zones for improved temperature control

Features

- Excellent mixing
- Wider range temperature control by heat fluid circulation
- Excellent heat exchange efficiency

Specification

Size LengthxWidthxThickness(mm)	Volume (μL)		Channel size Min. WidthxDepth(mm)
	Preheating	Mixing	
50x130x9.4	99	160	1x0.3



8 CYTOS[®]-200 series Residence time units

Laminar & Capillary RTU

Model numbers see below

Description

Two versions (laminar and capillary) are available for longer reaction times.

Reaction time is controlled by flow rate

Features

- Three different volumes
- Excellent heat exchange efficiency
- Optional Hastelloy versions are available for acidic liquids.
- Laminar type has plate structure.
- Capillary type has tubular structure.

Specification

Product name	Volume (mL)	Size LengthxWidthxThickness(mm)	Channel size Min. WidthxDepth(mm)	Material	Model number
Laminar0.9	0.9	50x130x8.2	2x0.2	SUS316	CYTOS200R0.9LCPC
Capillary0.9	0.9	50x130x30	φ0.6	SUS316	CYTOS200R0.9CPC
Capillary1.8	1.8	50x130x40	φ0.6	SUS316	CYTOS200R1.8CPC
Capillary4.5	4.5	50x130x55	φ0.9	SUS316	CYTOS200R4.5CPC
Laminar0.9	0.9	50x130x8.2	2x0.2	Hasttelloy	CYTOS200RHL0.9CPC
Capillary0.9	0.9	50x130x30	φ0.6	Hasttelloy	CYTOS200RH0.9CPC
Capillary1.8	1.8	50x130x40	φ0.6	Hasttelloy	CYTOS200RH1.8CPC
Capillary4.5	4.5	50x130x55	φ0.9	Hasttelloy	CYTOS200RH4.5CPC



Photo : Capillary1.8

Reactor (CYTOS[®]-L type)

Model numbers see below

Description


Original CPC design for CYTOS[®]-2000 reactor with internal heat fluid circulation zones for improved temperature control

Features

- Excellent heat exchange efficiency
- Optional Hastelloy versions available for acidic liquids.

Specification

Size LengthxWidthxThickness(mm)	Volume (μL)		Channel size Min. WidthxDepth(mm)	Material	Model number
	Preheating	Mixing			
100x150x10.8	740	1,200	1x0.3	SUS316	C2000-M-SUS
100x150x12.2	750	1,200	1x0.3	Hastelloy	C2000-M-HC

**Laminar & Capillary RTU**

Model numbers see below

Description

Two types (laminar and capillary) are available for longer reaction times

Features

- Two different volumes
- Excellent heat exchange efficiency
- Optional Hastelloy versions available for acidic liquids.
- Laminar type has plate structure.
- Capillary type has tubular structure.

Specification

Product name	Volume (mL)	Size LengthxWidthxThickness(mm)	Channel size Min. WidthxDepth(mm)	material	Model number
Laminar15	15	100x150x16.8	2x0.2	SUS316	CYTOS2000R15LCPC
Capillary15	15	100x150x25	φ2	SUS316	CYTOS2000R15CCPC
Capillary60	60	100x150x50	φ 2	SUS316	CYTOS2000R60CCPC
Laminar15	15	100x150x16.8	2x0.2	Hastelloy	CYTOS2000RHL15CCPC
Capillary15	15	100x150x25	φ 2	Hastelloy	CYTOS2000RH15CCPC
Capillary60	60	100x150x50	φ 2	Hastelloy	CYTOS2000RH60CCPC



Photo : Capillary15



11 KeyChem[®]-Lumino Reactor

Polaris (SUS)

Model number KCL-M-SUS

Description

Photoreaction reactor with 0.2 mm channel depth with excellent light irradiation efficiency

Features

- Excellent light irradiation efficiency
- Easy to be disassembled for cleaning

Specification

Size DiameterxThickness(mm)	Volume (μ L)	Channel size WidthxDepth(mm)
ϕ 37.5x0.5	110	1x0.2



12 KeyChem[®]-Lumino Reactors

Polaris (TiO₂)

Model numbers see below

Description

TiO₂ photocatalysis reactor

TiO₂ removal from the product is not necessary.

Features

- Suitable for photo TiO₂ catalysis experiments
- Can be disassembled for cleaning

Specification

Size DiameterxThickness(mm)	Volume (μ L)	Channel size WidthxDepth (mm)	Model number
ϕ 37.5x0.5	84	1x0.15	KCL-M-TiO ₂ -015
ϕ 37.5x0.5	170	1x0.3	KCL-M-TiO ₂ -030
ϕ 37.5x0.8	280	1x0.5	KCL-M-TiO ₂ -050





- Coupling type Thermo Stage
(Low temperature version / High temperature version)
- Stand alone type Thermo Stage
(Low temperature version / High temperature version)
- KeyChem[®]-Lumino Thermo Stage

C

Thermostatic units

1 KeyChem® series Thermostatic units

Model numbers see below

Coupling Type Thermo Stage (Low temperature version)

Description

Utilizes Peltier temperature control so no heat fluid is required


Several Thermo Stages can be combined for several reactions and residence time units

Features

- Peltier temperature control
- Temperature control range -15 - 80 °C
(for ambient temperature at 25 °C)
- Coupling of multiple units is possible.
- Oil-free operation (no heat fluid is not required.)

Specification

Dimension	(WxDxH) 270x230x220 mm(excluding projection)	Number of coupled Thermo Stages	Model number
Power requirements	AC 100 - 240 V, 50/60 Hz, 400 W (Maximum power requirement for 4 unit operation)	Coupling(2)	TS2-C
Temp. Range	Thermo Stage(M-L), (R-L) -15 - 80 °C (for ambient temperature at 25 °C)	Coupling(3)	TS3-C
Accuracy	±1 °C	Coupling(4)	TS4-C
Weight	Approx. 12 kg		



2 KeyChem® series Thermostatic units

Model numbers see below

Coupling Type Thermo Stage (High temperature version)

Description

Temperature control by electric heating

Heat fluid is not required.


Several Thermo Stages can be combined for several reactions and residence time units.

Features

- Temperature control by electric heating
- Temperature range, ambient temperature - 150 °C
- Coupling of multiple units is possible.
- Oil-free operation (no heat fluid is required.)

Specification

Dimension	(WxDxH) 270x230x255 mm (excluding projection)	Number of coupled Thermo Stages	Model number
Power requirements	AC 100 - 240 V, 50/60 Hz, 400 W (Maximum power requirement for 4 unit operation)	Coupling(2)	TS2-H
Temp. Range	Thermo Stage(M-H), (R-H) ambient temperature - 150 °C	Coupling(3)	TS3-H
Accuracy	±1 °C	Coupling(4)	TS4-H
Weight	Approx. 12 kg		





3 KeyChem[®] series Thermostatic units (for KeyChem[®]-Basic) Model numbers see below

Stand-alone Type Thermo Stage (Low temperature version)

Description

Peltier temperature control – no heat fluid required
 Several Thermo Stages can be coupled for multiple reactions
 Thermo Stage for RTU can also be added on.

Features

- Peltier temperature control
- Temperature control range -15 - 80 °C
 (for ambient temperature at 25 °C)
- Commercially available pumps can be connected.
- Oil-free operation (no heat fluid required.)



Specification

Dimension	(WxDxH) 120x230x220 mm (excluding projection)		Model number
Power requirements	AC 100 - 240 V, 50/60 Hz, 100 W	Mixer	TS-C-M
Type	Peltier temperature control	RTU	TS-C-R
Temp. Range	-15 - 80 °C (for ambient temperature at 25 °C)		
Accuracy	±1 °C		
Weight	Approx. 4 kg		



4 KeyChem[®] series thermostatic units (for KeyChem[®]-Basic) Model numbers see below

Stand-alone Type Thermo Stage (High temperature version)

Description

Temperature control by electric heater – no heat fluid required
 Several Thermo Stages can be coupled for multiple reactions
 Thermo Stage for RTU can also be added on.

Features

- Temperature control by electric heating
- Temperature range, ambient temperature - 150 °C
- Commercially available pumps can be connected.
- Oil-free operation (no heat fluid required.)



Specification

Dimension	(WxDxH) 120x230x255 mm (excluding projection)		Model number
Power requirements	AC 100 - 240 V, 50/60 Hz, 100 W	Mixer	TS-H-M
Temp. Range	Ambient temperature - 150 °C	RTU	TS-H-R
Accuracy	±1 °C		
Weight	Approx. 4 kg		



5 KeyChem[®] series Thermostat (for KeyChem[®]-lumino)
KeyChem-Lumino Thermo Stage

Model number TS-LU


Description

Dedicated thermostat for photoreaction.

A photoreaction reactor can be mounted

Features

- Peltier temperature control
- Any light sources can be used.
- Oil-free operation (no heat fluid required.)



Specification

Dimension	(WxDxH) 120x230x205 mm (excluding projection)
Power requirements	AC 100 - 240 V, 50/60 Hz, 100 W
Type	Peltier control
Temp. Range	-15 - 80 °C (for ambient temperature at 25 °C)
Accuracy	±1 °C
Weight	Approx. 4 kg



■ MR2 Syringe pump

D

Pump



MR2 Syringe pump

Model numbers see below


Description

Top quality syringe pump with all the functions that are required for microreaction

It can use several makes of syringes such as Hamilton, KLOEHN or disposable units. Since each syringe unit can be removed independently, any number of syringes can be set up for an experiment.

Features

- Easy to mount/remove different brands of syringes
- Any number of syringes can be set up for an experiment
- Continuous flow is obtained by alternating delivery from two syringes.
- Full PC control



Standard model

1-5 way valvex2

1-2 way valvex2

Specification

Dimension	(WxDxH) 230x350x280 (mm)
Power requirements	AC 100-240 V, 50/60 Hz, 300 W
Type	Syringe pump
Usable syringes	HAMILTON gas tight syringe(Luer lock type) Volume : 1, 2.5, 5, 10, 25 mL
	HSW disposable syringe(Luer lock type) Volume : 2(2.5), 5, 10, 20(18) mL *2, 20 mL The numbers in () are actual usage volumes.
	KLOEHN V6 syringe Volume : 2.5, 5, 10, 25 mL
Flow rate range	approx. 0.7 μ L/min - 55 mL/min *Depending on syringe size (recommended flow rate range : 10 μ L/min - 10 mL/min)
Flow rate accuracy	Within \pm 2%

Valve combination	Model number
1-5 way valvex4	MRSY04-40
1-5 way valvex2, 1-2 way valvex2	MRSY04-22



Syringe units

Can set up any number of syringes



With 1-5 way valvex2, 1-2 way valvex2

Rinsing can be done after reaction of two reagents.



With 1-5 way valvex4

Mixing of four reagents or continuous flow of two reagents is possible



- Connectors
- Ferrules
- Regulator
- Filters
- T type reactor
- Tubing
- Reactor housing
- Parts
- Others



Accessories - Consumable parts

1 Connectors

Tubing connector 1/16 inch (pack of 10) Stainless steel

Description

Stainless steel connector, to be used with ferrules (model numbers YMC-P-0011, YMC-P-0012) for connecting 1/16 inch tubing to a valve.

Material	Model number
SUS316	YMC-P-0001



2 Connectors

Tubing Connector 1/16 inch PEEK•M6HN

Description

PEEK connector to be used with ferrule (model number YMC-P-0013) for connecting 1/16 inch tubing to mixer housing.

Material	Model number
PEEK	YMC-P-0002



3 Connectors

Luer lock connector

Description

Luer lock connector for connecting a syringe to a valve

Material	Model number
PTFE / gland PCTFE / body	YMC-P-0003



4 Connectors

KeyChem-L Tubing connector (Low temperature version)

Description

Tubing connector for connecting 1/16 inch tubing to thermostat of KeyChem®-L (Low temperature version)

Material	Model number
POM	YMC-P-0004



5 Connectors

KeyChem[®]-L Tubing connector (High temperature version)

Description

Tubing connector for connecting 1/16 inch tubing to thermostat of KeyChem[®]-L (High temperature version)

Material	Model number
PEEK	YMC-P-0005



6 Connectors

Connector 1/16 inch (pack of 10) ETFE

Description

Connector for attaching 1/16 inch tubing to the syringe connectors (model number YMC-P-0009 and model number YMC-P-0010) with 1/16 inch tube.

Material	Model number
ETFE	YMC-P-0006



7 Connectors

EASYFIT10φ (pack of 10)

Description

Connector with built-in ferrule for 1/16 inch tubing.

Material	Model number
PEEK	YMC-P-0009



8 Connectors

Housing connector (Stainless steel tubing version)

Description

Connector for connecting 1/16 inch stainless steel tubing to the KeyChem[®] Mixer housing (Model number YMC-P-0030-01, and model number YMC-P-0030-02) requires ferrule (model number YMC-P-0001 and YMC-P-0011).

Material	Model number
SUS	YMC-P-0057



9 Adapter

Syringe adapter (female type)

Description

Adapter for connecting 1/16 inch tubing to a syringe (Luer lock type) using Connector 1/16 inch (Model number YMC-P-0006)

Material	Model number
ETFE	YMC-P-0007



10 Adapter

Luer adapter (male type)

Description

Adapter to be used with Connector 1/16 inch (Model number YMC-P-0006) for connecting a syringe needle for air sensitive reactions

Material	Model number
ETFE	YMC-P-0008



11 Ferrules

Ferrule 1/16 inch (pack par 10) PTFE

Description

Ferrule for use with connector (model number YMC-P-0001) for connecting tubing to a valve

Material	Model number
PTFE	YMC-P-0010



12 Ferrules

Ferrule 1/16 inch (pack par 10) Stainless steel

Description

Ferrule for use with connector (model number YMC-P-0001) for connecting tubing to a valve

Material	Model number
SUS316	YMC-P-0011



13 Ferrules

Ferrule 1/16 inch (pack par 10) PTFE•GF

Description

Ferrule for use with connector (model number YMC-P-0002, -0004, -0005, -0006) for connecting tubing to a valve

Material	Model number
PTFE/SUS316	YMC-P-0012



14-16 Back-pressure regulators

Back pressure regulator

Description

Regulator to be installed at the outlet of a reactor to increase pressure in the system

Material	Back pressure	Model number
PEEK	5 psi (0.03 MPa)	YMC-P-0013
	20 psi (0.14 MPa)	YMC-P-0014
	40 psi (0.28 MPa)	YMC-P-0015



17 In-line check valve

In-line check valve

Description

Check valve to be mounted at the inlet of a reactor to prevent reverse flow.

Material	Proof pressure	Model number
PEEK, gilding SUS, perfluorelastomer, tefzel	6.9 MPa	YMC-P-0052



18 Filters

In-line filter set (includes 10 filters)

Description

Filter to be mounted at the inlet of a reactor

Material	Size of bore	Model number
PTFE	25 µm	YMC-P-0016



19-20 Filters

Suction filter

Description

Suction filter for use on the suction side of a pump

Material	Model number
PP	YMC-P-0017
SUS	YMC-P-0018



21-22 For Microreactor

T type reactor (Stainless steel)

Description

T-tube reactor

I.D. (mm)	Material	Model number
0.5	SUS316	YMC-P-0019
1.0		YMC-P-0020



23-25 For Microreactor

T type reactor (PEEK)

Description

T-tube reactor made from PEEK

I.D. (mm)	Material	Model number
0.25	PEEK	YMC-P-0021
0.5		YMC-P-0022
1.0		YMC-P-0023



26-29 Connecting tubing

PTFE tubing

Description

Tubing for connecting various modules in a microreactor

I.D. (mm)	Material	O.D.	Length	Model number
0.25	PTFE	1/16 inch	10 m	YMC-P-0024
0.5				YMC-P-0025
0.75				YMC-P-0026
1.0				YMC-P-0027



30-31 Connecting tubing

Stainless steel tubing

Description

Tubing for connecting various modules in a microreactor

I.D. (mm)	Material	O.D.	Length	Model number
0.5	SUS316	1/16	10 m	YMC-P-0028
1.0		inch		YMC-P-0029



32 For KeyChem[®]

KeyChem[®] Mixer housing

Description

Mixer housing for KeyChem[®]

Simple experiments can be done with a KeyChem[®] mixer and commercially available thermostats such as oil or water baths.

One version is transparent for internal observation (YMC-P-0030-02)

Material	Model number
aluminum	YMC-P-0030-01
(A5052)	YMC-P-0030-02



YMC-P-0030-02

33 For KeyChem[®]

KeyChem[®] Mixer housing (Acrylic)

Description

Acrylic mixer housing for KeyChem[®] (A sister product of YMC-P-0030)

Transparent for internal observation

Material	Model number
Acrylic	YMC-P-0031



34 For KeyChem[®]

KeyChem[®] Residence time unit housing

Description

Housing for KeyChem[®] residence time units

Material	Model number
Aluminum (A5052)	YMC-P-0032



35 For CYTOS[®]-200 CYTOS[®]-200 Reactor housing

Description

Housing for CYTOS[®]-200 Reactor

Material	Model number
Aluminum (A5052)	YMC-P-0033



36 For CYTOS[®]-200 CYTOS[®]-200 Clamping unit

Description

Clamping unit for CYTOS[®]-200 Reactor

Model number
YMC-P-0034



37 For KeyChem[®]-Lumino KeyChem[®]-Lumino Reactor housing

Description

Housing for KeyChem[®]-Lumino Reactor

Part	Material	Model number
Cover	Aluminum (A5052), Quartz glass	YMC-P-0035
Gland	PTFE	
Base	SUS316	



38-41 For MR2 pump Hamilton Syringe attachment kit for MR2 pump

Description

Attachment kit for mounting a Hamilton syringe onto MR2 pump

Type	Model number
1 mL / 2.5 mL	YMC-P-0036
5 mL	YMC-P-0037
10 mL	YMC-P-0038
25 mL	YMC-P-0039



Photo : 25 mL type

42-45 For MR2 pump

Disposable Syringe attachment kit for MR2 pump

Description

Attachment kit for mounting a disposable syringe onto MR2 pump

Type	Model number
2 mL	YMC-P-0040
5 mL	YMC-P-0041
10 mL	YMC-P-0042
20 mL	YMC-P-0043



Photo : 20 mL

46 For MR2 pump

KLOEHN Syringe attachment kit for MR2 pump

Description

Attachment kit for mounting a KLOEHN syringe onto MR2 pump

Suitable for 2.5 mL, 5 mL, 10 mL, 25 mL syringes

Model number
YMC-P-0044




47 For MR2 pump

Reservoir base for MR2 pump

Description

Tray to be installed on the top of MR2 syringe pump for storing reagent bottles

Model number
YMC-P-0045



48 For MR2 pump

Reservoir base with pressure meter for MR2 pump

Description

Tray to be installed on the top of MR2 syringe pump for storing reagent bottles, complete with pressure meter

Model number
YMC-P-0046





49 For CYTOS®

O-ring for CYTOS® (Heating fluid line)

Description

O-ring for CYTOS® micro-reactor heat fluid tubing

Material	Model number
Fluoro-rubber (FPM)	YMC-P-0047



50 For CYTOS®

O-ring for CYTOS® (Reaction line)

Description

O-ring for CYTOS® micro-reactor reaction tubing

Material	Model number
Perfluoroelastomer	YMC-P-0048



51 Others

STACKING DRY BATH

Description

Oil-free dry heating bath for round bottomed flasks, suitable for 50, 100, 250, 500, 1000 mL flasks

Can be used on the standard 135 mm heater/stirrer

Conforms to JIS (Japanese Industrial Standards)

Model number
DRYBATH



52 For KeyChem[®]-Lumino UV light source unit (365 nm)

Model number YMC-P-0049

Description

Light Source for KeyChem[®]-Lumino with 6 LEDs in light emitting unit

Light type	Ultraviolet emitting LEDs
Peak wavelength	365 nm (25 °C)
Life time	500 hours in continuous operation
Light power	Max. 250 mW (each)
Power requirements	AC 100 - 240 V, 100 W
Dimensions	(WxDxH)
	Light emitting unit : 55x57x40 mm
	Power source unit : 50x250x161 mm



53 For KeyChem[®]-Lumino Low pressure mercury lamp for KeyChem[®]-Lumino(254 nm)

Model number YMC-P-0050

Description

Light source for KeyChem[®]-Lumino with 2 lamps in light emitting unit

Light type	Low pressure mercury lamp
Peak wavelength	254 nm
Life time	adv. 30000 hours
Light intensity	70 $\mu\text{W}/\text{cm}^2/10\text{ cm}$
Power requirements	AC 100 - 240 V, 18 W
Dimensions	(WxDxH)
	Light emitting unit : 57x100x45 mm
	Power source unit : 50x64x26.5 mm




54 For KeyChem[®]-Lumino Excimer lamp for KeyChem[®]-Lumino (308 nm)

Model number YMC-P-0051

Description

Excimer light source for KeyChem[®]-Lumino

Light type	Excimer lamp
Peak wavelength	308 nm
Life time	adv. 1000 hours
Power requirements	AC 100 V, 20 W
Dimensions	(WxDxH)
	Light emitting unit : 63x180x99 mm
	Power source unit : 50x250x161 mm



Please enquire about other lamp.