

Accessories for M-Series

Eyeieces M-series

- 10447159 Eyepiece 10x/21 adjustable
with eyegard
- 10447160 Eyepiece 10x/21B, adjustable
with eyegard
- 10445301 Eyepiece 16x/14B, adjustable
with eyegard
- 10445302 Eyepiece 25x/9.5B, adjustable
with eyegard
- 10445303 Eyepiece 40x/6B, adjustable
with eyegard

Measurement graticules for eyepieces


10376119 Graticule 12 mm : 120 divisions



10394771 Graticule 5 mm : 100 divisions




10376122 Graticule 100 x 1 sq. mm



10376120 Crosshair graticule



10447182 Graticule MZ16/MZ16 A/F/FA



10398408 Plain glass plate in mount

Stage micrometer

10310345 Stage micrometer, 50 mm
graduation marks 0,1 mm and 0,001 mm



Photo eyepieces

- 10446120 Photo eyepiece 8x
- 10445304 Photo eyepiece 10x
- 10445305 Photo eyepiece 16x
- 10445606 Spacing ring for photo eyepieces
(this space ring must be used in case of using a reticule for measurement
in the photo beam path together with the photo eyepiece)