

49000 - ET - DAPI

FILTER SET PARTS	PRICE
■ AT350/50x (EX)	\$150
■ T400lp (BS)	\$225
■ ET460/50m (EM)	\$325
Total Price	\$625.00

49021 - ET - EBFP2/Coumarin/Attenuated DAPI

Please see the excitation intensity note on page 6 for information on how to avoid rapid photobleaching of sample.

FILTER SET PARTS	PRICE
■ ET405/20x (EX)	\$325
■ T425lpxr (BS)	\$225
■ ET460/50m (EM)	\$325
Total Price	\$825.00

49028 - ET - DAPI for 395nm light sources

FILTER SET PARTS	PRICE
■ ET395/25x (EX)	\$325
■ T425lpxr (BS)	\$225
■ ET460/50m (EM)	\$325
Total Price	\$825.00

19011 - AT - Coumarin/Pacific Blue Longpass

FILTER SET PARTS	PRICE
■ AT405/30x (EX)	\$200
■ AT440DC (BS)	\$150
■ AT450lp (EM)	\$150
Total Price	\$450.00

TRF49902 - ET - 440nm Laser Bandpass Set for TIRF applications

For use with lasers between 430-445nm
Set includes additional longpass emission filter paired with bandpass emission filter in filter cube for greatly increased attenuation of TIRF lasers.

FILTER SET PARTS	PRICE
■ ZET436/20x (EX)	\$325
■ ZT442rdc (BS)	\$400
■ ET460lp (EM)	\$325
■ ET480/40m (EM)	\$325
Mounted in a Laser TIRF for Nikon TE2000/Ti	\$1,850
Mounted in a Laser TIRF for Olympus BX2 models	\$1,975
Mounted in a Laser TIRF for Zeiss Axio	\$1,725
Mounted in a Laser TIRF for Olympus BX3/IX3 models, for 25mm filters	\$1,975